

Dimensions:
Width: 6-3/4" (17,2 cm)
Length: 11-5/16" (28,7 cm)

REF EU-15955-CVT
Multi-Lumen Central Venous Catheterization Set
with Blue FlexTip® Catheter

ARROW®

Contents:

- 1:Multi-Lumen Indwelling Catheter: 9.5 Fr x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1:Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1:Introducer Needle: Echogenic 18 Ga x 2-1/2" (6.35 cm) XTW
- 1:Syringe: 5 ml Luer-Slip
- 1:Tissue Dilator 10 Fr (3.5 mm) x 4" (10.2 cm)
- 1:Fastener: Catheter Clamp
- 1:Clamp: Catheter
- 1:Scalpel: #11
- 1:ECG Extension Cable
- 5:Injection Site Cap

All components are CE 0086 unless otherwise noted.


Contains no medication.
Fluid path components are non-pyrogenic.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Многопросветный набор для центральной венозной катетеризации с катетерам Blue FlexTip
- HR** Višelumenski centralni venski kateter Set s Blue FlexTip katetrom
- CS** Vícelumenový centrální žilní katetizační set s katetrem Blue FlexTip
- DA** Multilumen centralvenøs kateteriserings-sæt med Blue FlexTip kateter
- NL** Meerlumen centraal veneuze katheterisatieset met Blue FlexTip katheter
- ET** Mitme valendikuga tsentraalse veenikateeteri komplekt koos Blue FlexTip kateetriga
- FI** Monilumeninen keskuskaskimokatetri sarja jossa Blue FlexTip katetri
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- EL** Σετ Καθητηριασμού Κεντρικής φλέβας Πολλαπλούς Αυλούς με Καθητήρα Blue FlexTip
- HU** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip katéterrel
- IS** Fjöllúmena holæðaleggsett með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- LT** Daugiakanalis centrinės venos kateterizavimo mini rinkinys su „Blue FlexTip“ kateteriu
- NO** Multilumen sentralt venekateriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanałowym z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- RO** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip
- RU** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- SR** Višelumenski centralni venski katetizacijski Set sa Blue FlexTip katetrom
- SK** Viacelúmenová centrálna venózná katetizačná súprava s katetrom Blue FlexTip
- SL** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- SV** Set för flerkanalig central venkateterisering med Blue FlexTip kateter
- TR** Multi Lumen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Багатопроевітний центральний венозний катетер Набор: катетер Blue с кінчиком Blue FlexTip Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (16 Ga)	0.48	3605
Medial 1 (14 Ga)	0.57	5390
Medial 2 (18 Ga)	0.38	1250
Medial 3 (18 Ga)	0.38	1160
Proximal (18 Ga)	0.41	1455

* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.20 ml.
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

5Z-159551-23A (5/14)

