

Dimensions:
Width: 6-3/4" (17,2 cm)
Length: 11-5/16" (28,7 cm)

REF EU-16402-N
Two-Lumen Central Venous Catheterization Set
with Blue FlexTip® Catheter
for Femoral Vein Insertion

ARROW®

Contents:

- 1: Two-Lumen Indwelling Catheter: 4 Fr x 12" (30 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: Catheter: 22 Ga x 1-3/4" (4.45 cm) Radiopaque over 25 Ga, RW Introducer Needle and 5 ml Luer-Slip Syringe
- 1: Introducer Needle: 20 Ga x 1-1/2" (4 cm) XTW
- 1: Tissue Dilator 5 Fr (1.8 mm) x 3" (7.5 cm)
- 1: Tissue Dilator 5 Fr (1.8 mm) x 1-5/8" (4.1 cm)
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Safety Scalpel: #11
- 1: ... SharpsAway II™ Locking Disposal Cup

All components are CE 0086 unless otherwise noted.


Contains no medication.
Fluid path components are non-pyrogenic.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REF Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Набор з дзьвухпросветным цэнтральным вярным катэтарам з наканечнікам Blue FlexTip для катэтарызацыі сцэгнавай вены
- HR** Komplet za centralnu vensku kateterizaciju s dva lumena i katetrom Blue FlexTip za uvođenje u femoralnu venu
- CS** Dvoulumenový centrální žilní katetizační set s katetrem Blue FlexTip, pro zavedení do femorální žíly
- DA** To-lumen centralvenekateriseringsæt med og Blue FlexTip kateter indføring i v. femoralis
- NL** Twee-lumen centraal veneuze katheterisatieset met Blue FlexTip catheter voor aanbrengen in de vena femoralis
- ET** Reieiveeni sisestamiseks mõeldud kahe valendikuga tsentraalveeni kateteriseerimise komplekt sinise FlexTip kateetriga
- FI** Kaksi-lumenoinen keskuslaskimokatetri sarja, jossa on sininen FlexTip katetri reisilaskimoon asettamiseen
- FR** Set de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip pour une insertion par la veine femorale
- DE** Zentralvenenkatheter-Set mit Doppel-Lumen mit Blue FlexTip Katheter zur Einführung in die V. femoralis
- EL** Σετ Καθητηριασμού Κεντρικής Φλέβας Διπλού Αυλού με Καθήτηρα Blue FlexTip για Εισαγωγή στη Μηριαία Φλέβα
- HU** Kétlumenes központi vénás katéterkészlet Blue FlexTip katéterrel a vena femoralis felőli bevezetéshez
- IS** Tveggja rása miðlæggt hölðaleggsett með Blue FlexTip æðalegg til ísetningar í lærleggsbláæð
- IT** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip per inserimento nella vena femorale
- LT** Dvieju kanalų centrinės venos kateterizavimo rinkinys su „Blue FlexTip“ kateteriu, skirtas įvesti į šlaunies veną
- NO** To-Lumen innsettingssett for centralt venekateter med Blue FlexTip kateter for innsetting i lårbensvene
- PL** Zestaw do cewnikowania żył centralnych z dwukanalowym cewnikiem z końcówką Blue FlexTip, do wprowadzenia przez żyłę udową
- PT** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip para Inserção na Veia Femoral
- RO** Set de cateterizare venoasă centrală cu două lumene cu cateter cu Blue FlexTip pentru introducerea în vena femorală
- RU** Набор с двухпросветным центральным венозным катетером с кончиком Blue FlexTip для катетеризации бедренной вены
- SR** Komplet sa dva lumena za centralnu vensku kateterizaciju sa plavim FlexTip katetrom za uvođenje u femoralnu venu
- SK** Dvojlúmenová centrálna venózna katetizačná súprava s katetrom Blue FlexTip na zavedenie do stehennej žily
- SL** Dvolumenski centralni venski pribor za kateterizacijo s katetrom Blue FlexTip za vstavljanje femoralno veno
- ES** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip para introducción en la vena femoral
- SV** Centralvenös katetersats med två lumen och blå FlexTip kateter för införing i vena femoralis
- TR** Femoral Ven Girişi için Blue FlexTip Kateter ile İki Lümenli Santral Venöz Kateterizasyon Seti
- UK** Набор з двохпросвітним центральним венозним катетером з кінчиком Blue FlexTip для катетеризації стегнової вени



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (22 Ga)	0.30	190
Proximal (22 Ga)	0.38	280

* Priming volumes are approximate and are done without the injection cap.
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

52-16402-100A (2/14)

LOT



(01) 0 08 01902 08775 1



REF EU-16402-N 4 Fr 2 Lumen 30 cm catheter length .018 inch dia. spring-wire guide
Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter for Femoral Vein Insertion

Peel to Open