

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

REF EU-16553-N
Multi-Lumen Central Venous Catheterization Set
with Blue FlexTip® Catheter

ARROW®

Contents:

- 1: Multi-Lumen Indwelling Catheter: 5.5 Fr x 5-1/8" (13 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1: Catheter: 22 Ga x 1-3/4" (4.45 cm) Radiopaque over 25 Ga, RW Introducer Needle and 5 ml Luer-Slip Syringe
- 1: Introducer Needle: 20 Ga x 1-1/2" (4 cm) XTW
- 1: Tissue Dilator 6 Fr (2.0 mm) x 1-5/8" (4.1 cm)
- 1: Tissue Dilator 6 Fr (2.0 mm) x 3" (7.6 cm)
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Safety Scalpel: #11
- 1: ... SharpsAway II™ Locking Disposal Cup

All components are CE 0086 unless otherwise noted.


Contains no medication.
Fluid path components are non-pyrogenic.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Многопросветный набор для центральной венозной катетеризации с катетером Blue FlexTip
- HR** Višelumenski centralni venski kateter Set s Blue FlexTip kateterom
- CS** Vícelumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** Multilumen centralvenøs kateteriseringssæt med Blue FlexTip kateter
- NL** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip katheter
- ET** Mitme valendikuga tsentraalse veenikateeteri komplekt koos Blue FlexTip kateetriga
- FI** Monilumeninen keskuslaskimokatetri sarja jossa sininen FlexTip katetri
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- EL** Σετ Καθητηριασμού Κεντρικής φλέβας Πολλαπλούς Αυλούς με Καθητήρα Blue FlexTip
- HU** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip katéterrel
- IS** Fjöllumena holæðaleggsett með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- LT** Daugiakanaliscentrinės venos kateterizavimo mini rinkinys su „Blue FlexTip“ kateteriu
- NO** Multilumen sentralt venekateriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanałowym z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- RO** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip
- RU** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- SR** Višelumenski centralni venski kateterizacioni Set sa Blue FlexTip kateterom
- SK** Viačlúmenová centrálna venózna katetrizačná súprava s katétrom Blue FlexTip
- SL** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- SV** Set för flerkanalig central venkateterisering med Blue FlexTip kateter
- TR** Multi Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Багатопрорітний центральний венозний катетер Набор: катетер Blue с кінчиком Blue FlexTip Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (20 Ga)	0.30	1200
Medial (22 Ga)	0.31	800
Proximal (22 Ga)	0.34	900

* Priming volumes are approximate and are done without the injection cap.
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

52-16553-100A (2/14)

LOT



(01) 0 08 01902 08776 8



REF EU-16553-N 5.5 Fr 3 Lumen 13 cm catheter length .018 inch dia. spring-wire guide
Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open