

# ARROW®



## Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CV-16702**    **7** Fr. | **2** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

### Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

**All components are CE 0086 unless otherwise noted.**

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**  
**Warning: Cancer and Reproductive Harm**  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Fluid path components are non-pyrogenic.



- fr** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip
- de** Besteck für zweilumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter
- it** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip
- pl** Dwukanałowy zestaw do cewnikowania żył centralnych z cewnikiem Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- ru** Комплект для катетеризации главных вен с двухпросветным катетером Blue FlexTip
- sl** Komplet za centralno vensko katetризacijo z dvema svetlinama s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip
- sv** Set för tvåkanalig central venkatetrering med Blue FlexTip-kateter
- tr** İki Lümenli Santral Venöz Kateterezasyon Seti, Blue FlexTip Kateter ile

| Catheter Cross Section | Lumen             | Priming Volume* (mL) | Flow Rate† (mL/hr) |
|------------------------|-------------------|----------------------|--------------------|
|                        |                   |                      |                    |
|                        | Distal (16 Ga.)   | 0.47                 | 3700               |
|                        | Proximal (16 Ga.) | 0.45                 | 3000               |

\* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

**EC REP** Teleflex Medical  
 IDA Business and Technology Park  
 Dublin Road, Athlone  
 Co. Westmeath, Ireland

**Manufacturer:**  
 Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

IV-16702-133C (01/18)



## Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CV-16702**    **7** Fr. | **2** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

# ARROW®