

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

REF **EU-16702-CVT**
Two-Lumen Central Venous Catheterization Set
with Blue FlexTip® Catheter
ARROW®

Contents:

- 1: Two-Lumen Indwelling Catheter:
7 Fr x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked:
.032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End — "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: Echogenic
18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Tissue Dilator: 8.5 Fr (2.8 mm) x 10.2 cm
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Scalpel: #11
- 1: ECG Extension Cable
- 1: Injection Site Cap

All components are CE 0086 unless otherwise noted.


Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Двухпросветный набор для центральной венозной катетеризации з катэтрам Blue FlexTip
- HR** Dvolumenski centralni venski kateter Set s Blue FlexTip kateterom
- CS** Dvoulumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** To-lumen centralvenøs kateteriseringsset med Blue FlexTip kateter
- NL** Centraal veneuze dubbellumen katheterisatieset met Blue FlexTip catheter
- ET** Kahe valendikuga tsentraalse veenikateetri komplekt koos Blue FlexTip kateetriga
- FI** Kaksiluumen-keskuslaskimokatetri setti jossa Blue FlexTip katetri
- FR** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip
- DE** Besteck für zweilumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας με Διπλούς Αυλούς και Καθετήρα με Blue FlexTip
- HU** Két lumenű centrális vénás katéterezési készlet Blue FlexTip katéterrel
- IS** Tveggja-lúmena holæðaleggsett með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip
- LT** Dviejų kanalų centrinės venos keteterizavimo rinkinys su „Blue FlexTip“ kateteriu
- NO** To-Lumen Sentralt venekateteriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żył centralnych z cewnikiem 2-kanalowym z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- RO** Trusa de cateterizare venoasă centrală, cu două lumene cu cateter Blue FlexTip
- RU** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip
- SR** Dvo-lumenski centralno venski set za kateterizaciju sa Blue FlexTip kateterom
- SK** Dvojlúmenová katetrizačná súprava na katetrizáciu centrálnej žily s katétrom Blue FlexTip
- SL** Dvolumenski centralni venski pribor za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip
- SV** Set för tvåkanalig central venkateterisering med Blue FlexTip kateter
- TR** Blue FlexTip Kateterli İki Lümenli Santral Venöz Kateterizasyon Seti
- UK** Двохпросвітний центральный венозный катетер з наконечником Blue FlexTip



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.47	3700
Proximal (16 Ga)	0.45	3000

*Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.
†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

5716702-1198 (10/14)

LOT



(01) 0 08 01902 06579 7



REF **EU-16702-CVT**
Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

7 Fr **2 Lumen** **20 cm** catheter length **.032** inch dia. spring-wire guide

Peel to Open