

Dimensions:  
Width: 6-3/4" (17.2 cm)  
Length: 11-5/16" (28.7 cm)

**REF EU-16702-EN**  
**Two-Lumen Central Venous Catheterization Set**  
with Blue FlexTip® Catheter  
**ARROW®**

**Contents:**

- 1: ..... Two-Lumen Indwelling Catheter: 7 Fr x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: ..... Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: ..... Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: ..... Arrow® Raulerson Spring-Wire Introduction Syringe: 5 ml
- 1: ..... Tissue Dilator: 8.5 Fr (2.8 mm) x 10.2 cm
- 1: ..... Fastener: Catheter Clamp
- 1: ..... Clamp: Catheter
- 1: ..... Pressure Transduction Probe

All components are CE 0086 unless otherwise noted.


Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

EC REP Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

- BE** Двухпросветный набор для центральной венозной катетеризации з катетерам Blue FlexTip
- HR** Dvolumenski centralni venski kateter Set s Blue FlexTip kateterom
- CS** Dvoulumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** To-lumen centralvenøs kateteriseringssæt med Blue FlexTip kateter
- NL** Centraal veneuze dubbellumen katheterisatieset met Blue FlexTip katheter
- ET** Kahe valendikuga tsentraalse veenikateetri komplekt koos Blue FlexTip kateetriga
- FI** Kaksiluumen-keskuslaskimokatetri setti jossa Blue FlexTip katetri
- FR** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip
- DE** Besteck für zweilumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας με Διπλούς Αυλούς και Καθετήρα με Blue FlexTip
- HU** Két lumenű centrális vénás katéterezési készlet Blue FlexTip katéterrel
- IS** Tveggja-lúmena holæðaleggsett með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip
- LT** Dviejų kanalų centrinės venos keteterizavimo rinkinys su „Blue FlexTip“ kateteriu
- NO** To-Lumen Sentralt venekateriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żył centralnych z cewnikiem 2-kanalowym z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- RO** Trusa de cateterizare venoasă centrală, cu două lumene cu cateter Blue FlexTip
- RU** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip
- SR** Dvo-lumenski centralno venski set za kateterizaciju sa Blue FlexTip kateterom
- SK** Dvojlúmenová katetrizačná súprava na katetrizáciu centrálnej žily s katétrom Blue FlexTip
- SL** Dvolumenski centralni venski pribor za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip
- SV** Set för tvåkanalig central venkatetrering med Blue FlexTip kateter
- TR** Blue FlexTip Kateterli İki Lümenli Santral Venöz Kateterizasyon Seti
- UK** Двохпросвітний центральний венозний катетер з наконечником Blue FlexTip



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.47	3700
Proximal (16 Ga)	0.45	3000

\*Priming volumes are approximate and are done without the injection cap.  
†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

57-16702-122B (11/14)

LOT



(01) 0 08 01902 07888 9



**REF EU-16702-EN** 7 Fr 2 Lumen 20 cm catheter length .032 inch dia. spring-wire guide  
**Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter**

Peel to Open