

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

REF EU-17702-CVT

Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

ARROW®

Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End — "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: Echogenic 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Tissue Dilator: 8.5 Fr (2.8 mm) x 10.2 cm
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Scalpel: #11
- 1: ECG Extension Cable
- 1: Injection Site Cap

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

BE Двухпросветный набор для центральной венозной катетеризации з катетэрам Blue FlexTip

HR Dvolumenski centralni venski kateter Set s Blue FlexTip kateterom

CS Dvolumenový centrální žilní katetrizační set s katetrem Blue FlexTip

DA To-lumen centralvenøs kateteriseringsæt med Blue FlexTip kateter

NL Centraal veneuze dubbellumen katheterisatieset met Blue FlexTip katheter

ET Kahe valendikuga tsentraalse veenikateetri komplekt koos Blue FlexTip kateetriga

FI Kaksiluumen-keskuslaskimokatetri setti jossa Blue FlexTip katetri

FR Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip

DE Besteck für zweilumigen, zentralen Venenkatheter mit Blue FlexTip Katheter

EL Σετ Καθετηριασμού Κεντρικής φλέβας με Διπλούς Αιχμούς και Καθετήρα με Blue FlexTip

HU Két lumenű centrális vénás katéterezési készlet Blue FlexTip katéterrel

IS Tveggja-lúmena holæðaleggsett með Blue FlexTip æðalegg

IT Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip

LT Dviejų kanalų centrinės venos keteterizavimo rinkinys su „Blue FlexTip“ kateteriu

NO To-Lumen Sentralt venekateteriseringssett med Blue FlexTip kateter

PL Zestaw do cewnikowania żył centralnych z cewnikiem 2-kanalowym z końcówką Blue FlexTip

PT Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip

RO Trusa de cateterizare venoasă centrală, cu două lumene cu cateter Blue FlexTip

RU Двухпросветный центральный венозный катетер с наконечником Blue FlexTip

SR Dvo-lumenski centralno venski set za kateterizaciju sa Blue FlexTip kateterom

SK Dvojlúmenová katetrizačná súprava na katetrizáciu centrálnej žily s katétrom Blue FlexTip


SL Dvolumenski centralni venski pribor za katetrizacijo s katetrom Blue FlexTip

ES Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip

SV Set för tvåkanalig central venkatetrisering med Blue FlexTip kateter

TR Blue FlexTip Kateterli İki Lümenli Santral Venöz Kateterizasyon Seti

UK Двохпросвітний центральный венозный катетер з наконечником Blue FlexTip



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (14 Ga)	0.57	5000
Proximal (18 Ga)	0.40	1500

*Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.


57-17702-1148 (10/14)

LOT



(01) 0 08 01902 06578 0

ARROW
INTERNATIONAL



REF

EU-17702-CVT

7^{Fr}

2 Lumen

20 cm catheter length

.032 inch dia. spring-wire guide

Peel to Open

Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter