

Dimensions:  
Width: 6-3/4" (17,2 cm)  
Length: 11-5/16" (28,7 cm)

**REF EU-17702-N**  
**Two-Lumen Central Venous Catheterization Set**  
with Blue FlexTip® Catheter  
**ARROW®**

**Contents:**

- 1: ..... Two-Lumen Indwelling Catheter:  
7 Fr x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: ..... Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: ..... Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: ..... Syringe: 5 ml Luer-Slip
- 1: ..... Tissue Dilator: 8.5 Fr (2.8 mm) x 10.2 cm
- 1: ..... Fastener: Catheter Clamp
- 1: ..... Clamp: Catheter

All components are CE 0086 unless otherwise noted.


Fluid path components are non-pyrogenic.  
Contains no medication.



Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

**EC REP** Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

- BE** Двухпросветный набор для центральной венозной катетеризации с катетерам Blue FlexTip
- HR** Dvolumenski centralni venski kateter Set s Blue FlexTip kateterom
- CS** Dvoulumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- DA** To-lumen centralvenøs kateteriseringsset med Blue FlexTip kateter
- NL** Centraal veneuze dubbellumen katheterisatieset met Blue FlexTip katheter
- ET** Kahe valendikuga tsentraalse veenikateetri komplekt koos Blue FlexTip kateetriga
- FI** Kaksiluumen-keskuslaskimokatetri setti jossa Blue FlexTip katetri
- FR** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip
- DE** Besteck für zweilumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας με Διπλούς Αυτούς και Καθετήρα με Blue FlexTip
- HU** Két lumenű centrális vénás katéterezési készlet Blue FlexTip katéterrel
- IS** Tveggja-lúmena holæðaleggsett með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip
- LT** Dviejų kanalų centrinės venos kateterizavimo rinkinys su „Blue FlexTip“ kateteriu
- NO** To-Lumen Sentralt venekateteriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żył centralnych z cewnikiem 2-kanalowym z końcówką Blue FlexTip
- PT** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- RO** Trusa de cateterizare venoasă centrală, cu două lumene cu cateter Blue FlexTip
- RU** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip
- SR** Dvo-lumenski centralno venski set za kateterizaciju sa Blue FlexTip kateterom
- SK** Dvojlúmenová katetrizačná súprava na katetrizáciu centrálnej žily s katétrom Blue FlexTip
- SL** Dvolumenski centralni venski pribor za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip
- SV** Set för tvåkanalig central venkateterisering med Blue FlexTip kateter
- TR** Blue FlexTip Kateterli İki Lümenli Santral Venöz Kateterizasyon Seti
- UK** Двохпросвітний центральний венозний катетер з наконечником Blue FlexTip



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (14 Ga)	0.57	5000
Proximal (18 Ga)	0.40	1500

\* Priming volumes are approximate and are done without the injection cap.  
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

SZ-17702-1248 (11/14)

**LOT**



(01) 0 08 01902 07889 6



**REF EU-17702-N**

**7 Fr**

**2 Lumen**

**20 cm**  
catheter length

**.032**  
inch dia. spring-wire guide

Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open