

Dimensions:
Width: 8-3/4" (22.2 cm)
Length: 11-1/8" (28.3 cm)

REF **EU-17702-CVTS**

Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

ARROW®

Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: Echogenic 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Tissue Dilator: 8.5 Fr (2.8 mm) x 10.2 cm
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Safety Scalpel: #11
- 1: SharpsAway II™ Locking Disposal Cup
- 1: ECG Extension Cable
- 1: Injection Site Cap

All components are CE 0086 unless otherwise noted.


Fluid path components are non-pyrogenic. Contains no medication.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Двухпросветный центральный венозный катетер з наконечником Blue FlexTip
- HR** Dvolumenski komplet za centralnu vensku kateterizaciju s plavim savitljivim vrhom FlexTip
- CS** Dvoulumenový centrální žilní katetizační set s katetrem Blue FlexTip
- DA** To-lumen centralvenekateteriseringssett med Blue FlexTip kateter
- NL** Centraal veneuze twee-lumen katheterisatieset met Blue FlexTip katheter
- ET** Kahe lumeniga tsentraalveeni kateteriseerimiskomplekt sinise painduva otsaga Blue kateetriga
- F** Kaksilumen-keskuslaskimokatetri setti, jossa Blue FlexTip katetri
- FR** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip
- DE** Besteck für zweilumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- EL** Σετ Καθετηριασμού Κεντρικής Φλέβας Διπλού Αυλού με Καθετήρα Blue FlexTip
- HU** Két lumenű centrális vénás katéterezési készlet Blue FlexTip katéterrel
- IS** Tveggja-lúmena sett fyrir holæðalegg með Blue FlexTip æðalegg
- IT** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip
- LT** Dviejų kanalų centrinės venos kateterizavimo rinkinys su „Blue FlexTip“ kateteriu
- NO** To-Lumen Sentralt venekateteriseringssett med Blue FlexTip kateter
- PL** Zestaw do cewnikowania żył centralnych, z dwukanałowym cewnikiem z końcówką Blue FlexTip i powłoką
- PT** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- RO** Trusa de cateterizare venoasă centrală cu două lumene, cu cateter Blue FlexTip
- RU** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip
- SR** Dvo-lumenski centralni venski kateterizacioni Set sa Blue FlexTip kateterom
- SK** Dvojlúmenová katetizačná súprava na katetrizáciu centrálnej žily s katetrom Blue FlexTip
- SL** Dvolumenski centralni venski pribor za katetrizacijo s katetrom Blue FlexTip
- ES** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip
- SV** Set för tvåkanalyg central venkateterisering med Blue FlexTip kateter
- TR** Blue FlexTip Kateterli İki Lümenli Santral Venöz Kateterizasyon Seti
- UK** Двохпросвітний центральний венозний катетер з наконечником Blue FlexTip



LOT

Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (14 Ga)	0.57	5000
Proximal (18 Ga)	0.40	1500

* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.


† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



(01) 0 08 01902 07353 2

SC17702:31A (12/14)

ARROW®
INTERNATIONAL



REF **EU-17702-CVTS**

7 Fr 2 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Peel to Open