



Multi-Lumen Antecubital Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CS-18763-E**

7 Fr. | **3** Lumen | **60** cm catheter length | **.035** inch dia. spring-wire guide

Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 60 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 39-1/2" (100 cm) PTFE-Coated (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 1: Dust Cap: Non-Vented
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Scalpel: #11
- 1: Stopcock: 3-way

All components are CE 0086 unless otherwise noted.

- fr** Jeu de cathétérisme veineux central antébrachial à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen, zentralen Venenkatheter zur Einführung in der Ellenbeuge mit Blue-FlexTip-Katheter
- it** Set per cateterismo della vena centrale cubitale anteriore multilume con catetere Blue FlexTip
- pl** Wielokanałowy zestaw do cewnikowania przedłokciowych żył centralnych z cewnikiem Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central Antecubital com Lúmen Múltiplo e Cateter Blue FlexTip
- ru** Комплект для катетеризации главных вен с многопросветным катетером Blue FlexTip для введения спереди от локтевого сустава
- sl** Komplet za antekubitalno centralno vensko katetrizacijo z več svetlinami s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central antecubital de luz múltiple con catéter Blue FlexTip
- sv** Set för flerkanalig antekubital central venkatetrisering med Blue FlexTip-kateter
- tr** Multi-Lümen Antekübital Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile



STERILE EO




Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65
Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Not made with natural rubber latex



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.70	1300
Medial (18 Ga.)	0.57	500
Proximal (18 Ga.)	0.57	500

* Priming volumes are approximate and are done without accessories.
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

DuPont™ Tyvek®
packaging

Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

CE REP
 0086
 Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

Teleflex®
www.teleflex.com

IS-18763-130C (01/18)



Multi-Lumen Antecubital Central Venous Catheterization Set with Blue FlexTip® Catheter

REF **CS-18763-E**

7 Fr. | **3** Lumen | **60** cm catheter length | **.035** inch dia. spring-wire guide