

Dimensions:
Width: 9-1/2" (24 cm)
Length: 12-1/2" (31.8 cm)

Contents:

- 1:Two-Lumen Indwelling Catheter: 7 Fr x 6" (16 cm) Pressure Injectable Radiopaque Polyurethane with Blue FlexTip[®], ARROWg^{ard}[®] Antimicrobial Surface Treatment[†], Extension Line Clamps
- 1:Spring-Wire Guide, Nitinol, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer[™]
- 1:Introducer Needle: Echogenic 18 Ga x 2-1/2" (6.35 cm) XTW
- 1:Arrow[®] Raulerson Spring-Wire Introduction Syringe: 5 ml
- 1:Tissue Dilator: 8.5 Fr (2.8 mm) x 10.2 cm
- 1:SecondSite[™] Adjustable Hub: Fastener
- 1:SecondSite[™] Adjustable Hub: Catheter Clamp
- 1:Safety Scalpel: #11
- 1:SharpsAway II[™] Locking Disposal Cup
- 1:Pressure Transduction Probe

[†] Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic. Not made with natural rubber latex.

Contraindications: The ARROWg^{ard} Blue[®] antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Store below 25°C (77°F).
Avoid excessive heat above 40°C (104°F).

Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA


EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

REF EU-22702-CVCPS

Pressure Injectable
**Two-Lumen Central Venous
Catheterization Kit, ARROWg^{ard} Blue[®]**

ARROW[®]

- BE** Набор «ARROWg^{ard} Blue» для центральной венозной катетеризации двух полостей органа, які ўпрыскаваецца пад націскам
- HR** Komplet za centralnu vensku kateterizaciju s dva lumena i tlačnim injektiranjem, ARROWg^{ard} Blue
- CS** Dvoulumenový centrální žilní katetrizační kit ARROWg^{ard} Blue, vhodný pro tlakové injekce
- DA** Tolumen centralvenekateteriseringsset til trykinjektion, ARROWg^{ard} Blue
- NL** Voor injectiepomp geschikte dubbel-lumen centraal veneuze katheterisatiekit, ARROWg^{ard} Blue
- ET** Survega sisestavut kahe valendikuga tsentraalveeni kateteriseerimise komplekt, ARROWg^{ard} Blue
- FI** Paineinjektoitava kaksilumeninen keskuskaskimokatetrisarja, ARROWg^{ard} Blue
- FR** Kit de cathétérisme veineux central à deux lumières, injectable sous pression, ARROWg^{ard} Blue
- DE** Zentralvenenkatheter-Kit mit Doppel-Lumen zur Druckinjektion, ARROWg^{ard} Blue
- EL** Kit Καθετηριασμού Κεντρικής Φλέβας Διπλού Αυλού με Δυνατότητα Έγχυσης με Πίεση, ARROWg^{ard} Blue
- HU** ARROWg^{ard} Blue nyomásalápoló befecskendezést alkalmazó kétlumenes központi vénás katéterkészlet
- IS** Tveggja rása miðlægt holleggjaset fyrir bláæðar sem er dælanlegt með þrýstingi, ARROWg^{ard} Blue
- IT** Kit per cateterismo a doppio lume per accesso venoso centrale iniettabile a pressione, ARROWg^{ard} Blue
- LT** Slėginio švirkštimo dviejų kanalų centrinės venos kateterizavimo rinkinys, „ARROWg^{ard} Blue“
- NO** Trykksatt injiserbart to-lumen sentralt venekateteriseringssett ARROWg^{ard} Blue
- PL** Dwukanałowy zestaw do cewnikowania żył centralnych ARROWg^{ard} Blue umożliwiający iniekcje pod wysokim ciśnieniem
- PT** Kit de Cateterização Venoso Central de Lúmen Duplo Injectável por Pressão ARROWg^{ard} Blue
- RO** Trusă de cateterizare venoasă centrală cu două lumene și injectare sub presiune, ARROWg^{ard} Blue
- RU** Двухпросветный набор для центральной венозной катетеризации с возможностью введения препарата под давлением, ARROWg^{ard} Blue
- SR** Komplet sa ubrizgavanjem pritiska i dva lumena za centralnu vensku kateterizaciju, ARROWg^{ard} Blue
- SK** Tlaková dvoulúmenová centrálna venózna katetrizačná súprava, ARROWg^{ard} Blue
- SL** Komplet za centralno vensko kateterizacijo in tlačno injiciranje z dvema kanaloma, ARROWg^{ard} Blue
- ES** Kit de cateterismo central venoso de dos lúmenes para inyección a presión ARROWg^{ard} Blue
- SV** Tryckinjicerbara centralvenösa katetersatser med två lumen, ARROWg^{ard} Blue
- TR** Basınçla Enjekte Edilebilen İki Lümenli Santral Venöz Kateterizasyon Seti, ARROWg^{ard} Blue
- UK** Двохпросвітний набір для центральної венозної катетеризації з можливістю введення препарату під тиском, ARROWg^{ard} Blue



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)	MAX Pressure Injection Flow Rate‡ (ml/sec)
Distal (18 Ga)	0.35	1300	5
Proximal (14 Ga)	0.50	5600	10

* Priming volumes are approximate and are done without injection caps.
† Based on 100 cm head height.
‡ Pressure injection flow rates are determined using media of 11.8 centipoise (cP) viscosity, with 60" administration set tubing.



SZ-22702-1090 (10/15)

REF EU-22702-CVCPS

7 Fr 2 Lumen 16 cm catheter length .025 inch dia. spring-wire guide

Pressure Injectable Two-Lumen Central Venous Catheterization Kit, ARROWg^{ard} Blue[®]

Peel to Open