

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

REF EU-22702-EN
Two-Lumen Central Venous Catheterization Set
with Blue FlexTip® ARROWg⁺ard Blue® Catheter

ARROW®

Contents:


- 1: Two-Lumen Indwelling Catheter:
7 Fr x 6" (16 cm) Radiopaque
Polyurethane with Blue FlexTip®,
ARROWg⁺ard® Antimicrobial Surface
Treatment[†], Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol,
Marked: .025" (0.64 mm) dia.
x 17-3/4" (45 cm) (Straight Soft Tip
on One End - "J" Tip on Other)
with Arrow Advancer™
- 1: Introducer Needle: 18 Ga x 2-1/2"
(6.35 cm) RW
- 1: Arrow® Raulerson Spring-Wire
Introduction Syringe: 5 ml
- 1: Tissue Dilator 8.5 Fr (2.8 mm)
x 4" (10.2 cm)
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Pressure Transduction Probe

[†] Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

Contraindications: The ARROWg⁺ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

Fluid path components are non-pyrogenic. Is not made with natural rubber latex.




Store below 25°C (77°F).
Avoid excessive heat
above 40°C (104°F).

Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REF Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Двухпросветный центральный венозный катетер з наконечником Blue FlexTip ARROWg⁺ard
- HR** Dvolumenski komplet za centralnu vensku kateterizaciju s plavim savitljivim vrhom FlexTip ARROWg⁺ard Blue
- CS** Dvoulumenový centrální žilní katetizační set s katetrem Blue FlexTip ARROWg⁺ard Blue
- DA** To-Lumen centralvenekateteriseringsssæt med Blue FlexTip ARROWg⁺ard Blue kateter
- NL** Centraal veneuze dubbellumen katheterisatieset met Blue FlexTip ARROWg⁺ard Blue katheter
- ET** Kahe luumeniga tsentraalveeni kateteriseerimiskomplekt sinise painduva otsaga ARROWg⁺ard Blue kateteriga
- FI** Kaksiluumen-keskuslaskimokatetri setti, jossa Blue FlexTip ARROWg⁺ard Blue katetri
- FR** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip ARROWg⁺ard Blue
- DE** Besteck für zweilumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg⁺ard Blue Katheter
- EL** Σετ Καθετηριασμού Κεντρικής Φλέβας Διπλού Αυλού με Καθετήρα Blue FlexTip ARROWg⁺ard Blue
- HU** Két lumenű centrális vénás katéterezési készlet Blue FlexTip ARROWg⁺ard Blue katéterrel
- IS** Tveggja-lúmena sett fyrir holæðalegg með Blue FlexTip ARROWg⁺ard Blue æðalegg
- IT** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip ARROWg⁺ard Blue
- LT** Dviejų kanalų centrinės venos kateterizavimo rinkinys su „Blue FlexTip ARROWg⁺ard Blue“ kateteriu
- NO** To-Lumen Sentralt venekateteriseringssett med Blue FlexTip ARROWg⁺ard Blue kateter
- PL** Zestaw do cewnikowania żył centralnych, z 2-kanalowym cewnikiem z końcówką Blue FlexTip i powłoką ARROWg⁺ard Blue
- PT** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip, ARROWg⁺ard Blue
- RO** Trusa de cateterizare venoasă centrală cu două lumene, cu cateter Blue FlexTip ARROWg⁺ard Blue
- RU** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip ARROWg⁺ard
- SR** Dvo-lumenski centralni venski kateterizacioni Set sa Blue FlexTip ARROWg⁺ard Blue kateterom
- SK** Dvojlúmenová katetizačná súprava na katetrizáciu centrálnej žily s katetrom Blue FlexTip ARROWg⁺ard Blue
- SL** Dvolumenski centralni venski pribor za katetrizacijo s katetrom Blue FlexTip ARROWg⁺ard Blue
- ES** Conjunto de cateterización venosa central de dos lúmenes con catéter Blue FlexTip ARROWg⁺ard Blue
- SV** Set för tvåkanalyg central venkateterisering med Blue FlexTip ARROWg⁺ard Blue kateter
- TR** Blue FlexTip ARROWg⁺ard Blue Kateterli İki Lümenli Santral Venöz Kateterizasyon Seti
- UK** Двохпросвітний центральный венозный катетер з наконечником Blue FlexTip ARROWg⁺ard



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (18 Ga)	0.32	1500
Proximal (14 Ga)	0.47	5000

*Priming volumes are approximate and are done without the injection cap.
†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

SI22702-111A (2/14)

LOT



(01) 0 08 01902 08819 2



REF EU-22702-EN 7 Fr 2 Lumen 16 cm catheter length .025 inch dia. spring-wire guide
Two-Lumen Central Venous Catheterization Set with ARROWg⁺ard Blue® Catheter

Peel to Open