



REF CV-22703-E

ARROW

7 Fr. | 3 Lumen | 16 cm catheter length | .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg⁺ard Blue® Catheter

Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg⁺ard® Antimicrobial Surface Treatment¹, Extension Line Clamps
 - 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer
 - 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
 - 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
 - 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
 - 1: SecondSite™ Adjustable Hub: Fastener
 - 1: SecondSite™ Adjustable Hub: Catheter Clamp
 - 3: Dust Cap: Non-Vented
 - 1: Pressure Transduction Probe
- ¹ Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

⚠ California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.


Not made with natural rubber latex.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The ARROWg⁺ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip ARROWg⁺ard Blue
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg⁺ard Blue Katheter
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip ARROWg⁺ard Blue
- pl** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanałowym z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg⁺ard Blue
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip, ARROWg⁺ard Blue
- ru** Набор с многопросветным центральным венозным катетером ARROWg⁺ard Blue с кончиком Blue FlexTip
- sl** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip ARROWg⁺ard Blue
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip ARROWg⁺ard Blue
- sv** Set för flerkanalig central venkatetrisering med Blue FlexTip ARROWg⁺ard Blue kateter
- tr** Multi Lumen Santral Venöz Kateterizasyon Seti, Mavi esnek uçlu ARROWg⁺ard Blue Kateterli



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.39	3400
Medial (18 Ga.)	0.35	1800
Proximal (18 Ga.)	0.37	1900

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



LOT

Use
ByEC REP
CE
0086

Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

SZ-22703-104I (10/17)



REF CV-22703-E

7 Fr. | 3 Lumen | 16 cm catheter length | .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip®
ARROWg⁺ard Blue® Catheter

ARROW