

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

REF EU-22703-EN
Multi-Lumen Central Venous Catheterization Set
with Blue FlexTip® ARROWg⁺ard Blue® Catheter

ARROW®

Contents:

- 1:Multi-Lumen Indwelling Catheter:
7 Fr x 6" (16 cm) Radiopaque Polyurethane with Blue FlexTip®, ARROWg⁺ard® Antimicrobial Surface Treatment[†], Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1:Arrow® Raulerson Spring-Wire Introduction Syringe: 5 ml
- 1: Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW
- 1:Tissue Dilator 8.5 Fr (2.8 mm) x 4" (10.2 cm)
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Pressure Transduction Probe

[†] Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

Contraindications: The ARROWg⁺ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

Fluid path components are non-pyrogenic. Is not made with natural rubber latex.




Store below 25°C (77°F).
Avoid excessive heat above 40°C (104°F).

Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Многопросветный центральный венозный катетер Набор: Катетер Blue с наконечником Blue FlexTip ARROWg⁺ard Blue Catheter
- HR** Komplet za višelumensku centralno vensku kateterizaciju sa Blue FlexTip ARROWg⁺ard Blue kateterom
- CS** Všelumenový centrální žilný katetizační set s katetrem Blue FlexTip ARROWg⁺ard Blue
- DA** Multilumen centralvenøs kateteriseringsæt med Blue FlexTip ARROWg⁺ard Blue catheter
- NL** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip ARROWg⁺ard Blue katheter
- ET** Mitme valendikuga tsentraalse veenikateetri komplekt koos Blue FlexTip ARROWg⁺ard Blue kateetriga
- FI** Moniluumeninen keskuskaskimokatetri sarja, jossa sininen FlexTip ARROWg⁺ard Blue katetri
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip ARROWg⁺ard Blue
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg⁺ard Blue Katheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας με Πολλαπλούς Αιρούς και Καθετήρα ARROWg⁺ard Blue με Blue FlexTip
- HU** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip ARROWg⁺ard Blue katéterrel
- IS** Fjöllúmena holæðaleggsett með Blue FlexTip ARROWg⁺ard Blue æðalegg
- IT** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip ARROWg⁺ard Blue
- LT** Daugiakanaliscentrinės venos kateterizavimo mini rinkinys su „Blue FlexTip ARROWg⁺ard Blue“ antimikrobinio kateteriu
- NO** Multilumen Sentralt venekateriseringssett med Blue FlexTip ARROWg⁺ard Blue kateter
- PL** Zestaw do cewnikowania żyły centralnej z cewnikami wielokanalowym z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg⁺ard Blue
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip, ARROWg⁺ard Blue
- RO** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip ARROWg⁺ard Blue
- RU** Набор с многопросветным центральным венозным катетером ARROWg⁺ard Blue с кончиком Blue FlexTip
- SR** Višelumenski centralni venski kateterizacioni Set sa Blue FlexTip ARROWg⁺ard Blue kateterom
- SK** Viacľúmenová centrálna venózna katetizačná súprava s katetrom Blue FlexTip ARROWg⁺ard Blue
- SL** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip ARROWg⁺ard Blue
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip ARROWg⁺ard Blue
- SV** Set för flerkanalig central venkatetrering med Blue FlexTip ARROWg⁺ard Blue kateter
- TR** ARROWg⁺ard Blue Multi Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Багатопроектний центральний венозний катетер Набор: катетер Blue с кінчиком Blue FlexTip ARROWg⁺ard Blue Catheter



| Lumen | Priming Volume* (ml) | Flow Rate† (ml/hour) |
|------------------|----------------------|----------------------|
| Distal (16 Ga) | 0.39 | 3400 |
| Medial (18 Ga) | 0.35 | 1800 |
| Proximal (18 Ga) | 0.37 | 1900 |

* Priming volumes are approximate and are done without the injection cap.
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

SE22703-16A 02/14

LOT



(01) 0 08 01902 08820 8



REF EU-22703-EN 7 Fr 3 Lumen 16 cm catheter length .032 inch dia. spring-wire guide
Multi-Lumen Central Venous Catheterization Set with ARROWg⁺ard Blue® Catheter

Peel to Open