



REF CV-22854



8.5 Fr. | 4 Lumen | 16 cm catheter length | .032 inch dia. spring-wire guide

# Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg<sup>+</sup>ard Blue® Catheter

### Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr. x 16 cm radiopaque Polyurethane with Blue FlexTip®, ARROWg<sup>+</sup>ard® Antimicrobial Surface Treatment<sup>1</sup>, Extension Line Clamps
  - 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
  - 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
  - 1: Syringe: 5 mL Luer-Slip
  - 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
  - 1: SecondSite™ Adjustable Hub: Fastener
  - 1: SecondSite™ Adjustable Hub: Catheter Clamp
  - 4: Dust Cap: Non-Vented
- <sup>1</sup> Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**

**Warning: Cancer and Reproductive Harm** [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).


Not made with natural rubber latex.

Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

**Contraindications: The ARROWg<sup>+</sup>ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).**

- fr** Jeu de cathétérisme veineux central à quatre lumières avec cathéter Blue FlexTip ARROWg<sup>+</sup>ard Blue
- de** Besteck für vierlumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg<sup>+</sup>ard Blue Katheter
- it** Set per cateterismo della vena centrale a quattro lumi con catetere Blue FlexTip ARROWg<sup>+</sup>ard Blue
- pl** Zestaw do cewnikowania żyły centralnej z 4-kanalowym cewnikiem z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg<sup>+</sup>ard Blue
- pt** Conjunto de Cateterização Venosa Central com Lúmen Quádruplo e Cateter Blue FlexTip, ARROWg<sup>+</sup>ard Blue
- ru** Четырех-просветный набор для центральной венозной катетеризации с катетером Blue FlexTip ARROWg<sup>+</sup>ard Blue
- sl** Štiri-lumenski centralni venski komplet za katetризacijo s katetrom Blue FlexTip ARROWg<sup>+</sup>ard Blue
- es** Conjunto de cateterización venosa central de cuatro lúces con catéter Blue FlexTip ARROWg<sup>+</sup>ard Blue
- sv** Set för fyrkanalig central venkatetrisering med Blue FlexTip ARROWg<sup>+</sup>ard Blue kateter
- tr** Dört Lümenli Santral Venöz Kateterizasyon Seti, Mavi esnek uçlu ARROWg<sup>+</sup>ard Blue Kateterli



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.39	3200
Medial 1 (14 Ga.)	0.47	5400
Medial 2 (18 Ga.)	0.36	1600
Proximal (18 Ga.)	0.38	1800

\* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



STERILE EO



LOT



CE REP 0086

Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Manufacturer:  
Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

SZ-22854-1001 (10/17)



REF CV-22854

8.5 Fr. | 4 Lumen | 16 cm catheter length | .032 inch dia. spring-wire guide

Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg<sup>+</sup>ard Blue® Catheter

