

Dimensions:
 Width: 6-3/4" (17,2 cm)
 Length: 11-5/16" (28,7 cm)

Antimicrobial surfaces utilizing
 chlorhexidine and silver sulfadiazine ‡

REF BR-22854

Quad-Lumen Central Venous Catheterization Set

with Blue FlexTip® ARROWg^{ard} Blue® Catheter

ARROW®

Special Order Product

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr x 6" (16 cm) Radiopaque Polyurethane with Blue FlexTip®, ARROWg^{ard}® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™
- 1: Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Catheter: 18 Ga x 2-1/2" (6.35 cm) Radiopaque over 20 Ga RW Introducer Needle with 5 ml Luer-Slip Syringe
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Scalpel: #11
- 4: Stopcocks: 4-way High-Flow
- 1: Tissue Dilator

¹ Licensed under US patent No. 6,706,024.

‡ **Contraindications: The ARROWg^{ard} Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).**


Fluid path components are non-pyrogenic.
 Does not contain natural rubber latex.



Store below 25°C (77°F).
 Avoid excessive heat
 above 40°C (104°F).

Arrow International, Inc.
 2400 Bernville Road
 Reading, PA 19605 USA

EC REP Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (16 Ga)	0.39	3200
Medial 1 (14 Ga)	0.47	5400
Medial 2 (18 Ga)	0.36	1600
Proximal (18 Ga)	0.38	1800

* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

LOT



(01) 0 08 01902 03759 6

SZ:22854-104E (11/12)

ARROW
 INTERNATIONAL
 Special Order Product

REF BR-22854 8.5 Fr 4 Lumen 16 cm catheter length .032 inch dia. spring-wire guide
 Quad-Lumen Central Venous Catheterization Set with ARROWg^{ard} Blue® Catheter

Peel to Open