

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

REF EU-22854-N

Quad-Lumen Central Venous Catheterization Set

with Blue FlexTip® ARROWg⁺ard Blue® Catheter

ARROW®

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr x 6" (16 cm) Radiopaque Polyurethane with Blue FlexTip®, ARROWg⁺ard® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™
- 1: Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Tissue Dilator: 9 Fr (3.0 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

¹ Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic. Not made with natural rubber latex.

Contraindications: The ARROWg⁺ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Store below 25°C (77°F).
Avoid excessive heat above 40°C (104°F).

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Teleflex®
www.teleflex.com


LOT

- BE** Чатырохпросветны цэнтральны венозны катэтер Набор: Катэтер Blue с наконечником Blue FlexTip ARROWg⁺ard Blue Catheter
- HR** Quad-Lumen Komplet za centralnu vensku kateterizaciju s Blue FlexTip ARROWg⁺ard Blue kateterom
- CS** Čtyřlumenový centrální žilní katetizační set s katetrem Blue FlexTip ARROWg⁺ard Blue
- DA** Firdobbelt lumen centralvenekateteriseringsset med Blue FlexTip ARROWg⁺ard Blue kateter
- NL** Centraal veneuze quad-lumen katheterisatieset met Blue FlexTip ARROWg⁺ard Blue katheter
- ET** Neljavalendikulisea tsentraalse veenikateeteri komplekt koos Blue FlexTip ARROWg⁺ard Blue kateetriga
- FI** Quad-Lumen keskuslaskimokatetri sarja jossa sininen FlexTip ARROWg⁺ard Blue katetri
- FR** Jeu de cathétérisme veineux central à quatre lumières avec cathéter Blue FlexTip ARROWg⁺ard Blue
- DE** Besteck für vierlumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg⁺ard Blue Katheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας με τετραπλάσιο αυλού και Καθετήρα ARROWg⁺ard Blue με Blue FlexTip
- HU** Quad-Lumen centrális vénás katéterezési szett, Blue FlexTip ARROWg⁺ard Blue katéterrel
- IS** Fjögur lúmena holæðaleggsett með Blue FlexTip ARROWg⁺ard Blue æðalegg
- IT** Set per caterismo della vena centrale a quattro lumi con catetere Blue FlexTip ARROWg⁺ard Blue
- LT** Keturių angų centrinės venos kateterizavimo mini rinkinys su „Blue FlexTip ARROWg⁺ard Blue“ antimikrobinium kateterium
- NO** Sentralt venekateteriseringssett med fire lumen og Blue FlexTip ARROWg⁺ard Blue-kateter
- PL** Zestaw do cewnikowania żyły centralnej z 4-kanalowym cewnikiem z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg⁺ard Blue
- PT** Conjunto de Cateterização Venosa Central com Lúmen Quádruplo e Cateter Blue FlexTip, ARROWg⁺ard Blue
- RO** Set de cateterizare venoasă centrală quad-lumen cu cateter Blue FlexTip ARROWg⁺ard Blue
- RU** Четырех-просветный набор для центральной венозной катетеризации с катетером Blue FlexTip ARROWg⁺ard Blue
- SR** Quad-Lumen centralni venski kateterizacioni Set sa Blue FlexTip ARROWg⁺ard Blue kateterom
- SK** Štvorlúmenová centrálna venózna katetizačná súprava s katétrom Blue FlexTip ARROWg⁺ard Blue
- SL** Štiri-lumenski centralni venski komplet za katetizacijo s katetrom Blue FlexTip ARROWg⁺ard Blue
- ES** Conjunto de cateterización venosa central de cuatro luces con catéter Blue FlexTip ARROWg⁺ard Blue
- SV** Set för fyrkanalyg central venkatetrering med Blue FlexTip ARROWg⁺ard Blue kateter
- TR** ARROWg⁺ard Blue Dörtlü-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Чотирипросвітний центральний венозний катетер Набор: катетер Blue с кінчиком Blue FlexTip ARROWg⁺ard Blue Catheter

Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.39	3200
Medial 1 (14 Ga)	0.47	5400
Medial 2 (18 Ga)	0.36	1600
Proximal (18 Ga)	0.38	1800

*Priming volumes are approximate and are done without the injection cap.

†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Catheter Cross Section

5Z22854-111A (1/16)

ARROW®
INTERNATIONAL



REF EU-22854-N

8.5 Fr 4 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Quad-Lumen Central Venous Catheterization Set with ARROWg⁺ard Blue® Catheter

Point to Order