

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

REF EU-22855-EN

Multi-Lumen Central Venous Catheterization Set

with Blue FlexTip® ARROWg⁺ard Blue® Catheter

ARROW®

Contents:

- 1: Multi-Lumen Indwelling Catheter: 8.5 Fr x 6" (16 cm) Radiopaque Polyurethane with Blue FlexTip®, ARROWg⁺ard Antimicrobial Surface Treatment[†], Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1: Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 ml
- 1: Tissue Dilator 9 Fr (3.0 mm) x 4" (10.2 cm)
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Pressure Transduction Probe

[†] Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

Contraindications: The ARROWg⁺ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

Is not made with natural rubber latex. Fluid path components are non-pyrogenic.




Store below 25°C (77°F).
Avoid excessive heat above 40°C (104°F).

Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Многопросветный центральный венозный катетер Набор: Катетер Blue с наконечником Blue FlexTip ARROWg⁺ard Blue Catheter
- HR** Komplet za višelumensku centralno vensku kateterizaciju sa Blue FlexTip ARROWg⁺ard Blue kateterom
- CS** Vícelumenový centrální žilní katetizační set s katetrem Blue FlexTip ARROWg⁺ard Blue
- DA** Multilumen centralvenøs kateteriseringssæt med Blue FlexTip ARROWg⁺ard Blue kateter
- NL** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip ARROWg⁺ard Blue katheter
- ET** Mitme valendikuga tsentraalse veenikateetri komplekt koos Blue FlexTip ARROWg⁺ard Blue kateetriga
- FI** Monilumenninen keskuslaskimokatetri sarja, jossa sininen FlexTip ARROWg⁺ard Blue katetri
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip ARROWg⁺ard Blue
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg⁺ard Blue Katheter
- EL** Σετ Καθητηριασμού Κεντρικής Φλέβας με Πολλαπλούς Αιλούς και Καθέτηρα ARROWg⁺ard Blue με Blue FlexTip
- HU** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip ARROWg⁺ard Blue katéterrel
- IS** Fjöllumena holæðaleggsett með Blue FlexTip ARROWg⁺ard Blue æðalegg
- IT** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip ARROWg⁺ard Blue
- LT** Daugiakanalis centrinės venos kateterizavimo mini rinkinys su „Blue FlexTip ARROWg⁺ard Blue“ antimikrobinu kateteriu
- NO** Multilumen Sentralt venekateriseringssett med Blue FlexTip ARROWg⁺ard Blue kateter
- PL** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanałowym z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg⁺ard Blue
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip, ARROWg⁺ard Blue
- RO** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip ARROWg⁺ard Blue
- RU** Набор с многопросветным центральным венозным катетером ARROWg⁺ard Blue с кончиком Blue FlexTip
- SR** Višelumenski centralni venski kateterizacijski set sa Blue FlexTip ARROWg⁺ard Blue kateterom
- SK** Viačlúmenová centrálna venózna katetizačná súprava s katétrom Blue FlexTip ARROWg⁺ard Blue
- SL** Večlumenski centralni venski komplet za katetizacijo s katetrom Blue FlexTip ARROWg⁺ard Blue
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip ARROWg⁺ard Blue
- SV** Set för flerkanalig centralvenkateterisering med Blue FlexTip ARROWg⁺ard Blue kateter
- TR** ARROWg⁺ard Blue Multi Lumen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Багатопрорітний центральний венозний катетер Набор: катетер Blue FlexTip ARROWg⁺ard Blue Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (16 Ga)	0.58	3650
Medial 1 (14 Ga)	0.71	6590
Medial 2 (18 Ga)	0.36	1260
Medial 3 (18 Ga)	0.38	1430
Proximal (18 Ga)	0.40	1460

*Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

5222855-107A (4/14)



(01) 0 08 01902 08858 1



REF EU-22855-EN 8.5 Fr 5 Lumen 16 cm catheter length .032 inch dia. spring-wire guide
Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg⁺ard Blue® Catheter

Peel to Open