

Dimensions:
Width: 10-5/8" (27 cm)
Length: 12-1/2" (31.8 cm)

Contents:

- 1:Indwelling Catheter: 4 Fr x 15-3/4" (40 cm)
Pressure Injectable Radiopaque Polyurethane
with Blue FlexTip®, T-Port Connector, Extension Line
Clamps, Contamination Guard and Placement Wire
- 1:Spring-Wire Guide, Nitinol, Marked: .018"
(0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft
Tip on One End - Straight Stiff Tip on Other)
with Arrow Advancer™
- 1:Introducer Needle: Echogenic 21 Ga x 2-3/4"
(7.00 cm) TW
- 1:Syringe: 10 ml Luer-Lock
- 1:Fastener: Catheter Clamp
- 1:Clamp: Catheter
- 1:GlideThru™ Peel-Away Sheath: 4 Fr x 2-3/4"
(7 cm) Radiopaque over 4 Fr Dilator
- 1:CSR Wrap
- 1:Dressing: GRIP-LOK® Catheter Fixation Device
- 1:Safety Scalpel: #11
- 1:Catheter Trimmer
- 1:SharpsAway II™ Locking Disposal Cup
- 1:Tourniquet²
- 1:Paper Tape Measure²
- 1:Valve: MicroCLAVE® Neutral Displacement Connector
- 1:Flow Rate Information Card with Injection Log²
- 1:Patient ID Card²
- 1:Chart Sticker²
- 1:Patient Information Booklet²

¹ A registered trademark of ICU Medical, Inc.

² Packaged externally to the sterile kit.

³ A registered trademark of TIDI

Fluid path components are non-pyrogenic.
Contains no medication.

All components are CE 0086 unless otherwise noted.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

REF EU-24041-HP

**Pressure Injectable
PICC Kit
with Blue FlexTip® Catheter**

ARROW®

NL Kit met PICC voor hogedrukinjectie,
met Blue FlexTip catheter

FR Kit PICC pour injection sous pression
avec cathéter Blue FlexTip

DE Pressure Injectable PICC (Peripherally
Inserted Central Catheter) Kit
mit Blue FlexTip Katheter

EL Κιτ PICC για έγχυση υπό πίεση,
που περιλαμβάνει καθετήρα με άκρο
Blue FlexTip


IT Kit con PICC per infusione a pressione
con catetere Blue FlexTip

PL Zestaw nadającego się do wstrzykiwania
pod ciśnieniem cewnika PICC z końcówką
Blue FlexTip

PT Kit de cateter central de inserção periférica
para injeção a pressão com cateter
com ponta Blue FlexTip

ES Kit de cateterización central de inserción
periférica para inyección a presión
con catéter con punta Blue FlexTip

SV Set med tryckinjicerbar perifert införd
central kateter med Blue FlexTip



LOT

Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (18 Ga)	0.40	1020	4679	4

* Priming volumes are approximate and are done without injection caps.

† Flow rates were determined using room temperature water, 100 cm head height

and represent approximate flow capabilities.

‡ Pump flow rates are determined at maximum pump pressure of 10 psig and represent

approximate flow capabilities.

**Pressure injection flow rates are determined at the injector pressure setting of 300 psi

maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.



(01) 0 08 01902 09289 2

52-24041-103A (3/15)

ARROW®
INTERNATIONAL



REF

EU-24041-HP

Pressure Injectable PICC Kit with Blue FlexTip® Catheter

4 Fr

40 cm

catheter
length

.018 inch dia.
spring-wire
guide

Radio
opaque