

Dimensions:
Width: 11" (28 cm)
Length: 19" (48.3 cm)




Arrow® Pressure Injectable Two-Lumen PICC Kit

REF EU-24052-HPMSB

Contents:

- | | |
|--|--|
| <ul style="list-style-type: none"> 1: Two-Lumen Indwelling Catheter: 5 Fr x 15-3/4" (40 cm) Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, T-Port Connector, Extension Line Clamps, Contamination Guard and Placement Wire 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer™ 1: Filter Needle: 18 Ga TW x 1-1/2" (3.81 cm) with Blunt Fill Tip 1: Introducer Needle: Echogenic 21 Ga x 2-3/4" (7.00 cm) IW 1: Eclipse™ Safety Needle: 22 Ga x 1-5/8" (4 cm) 1: Eclipse™ Safety Needle: 25 Ga x 1" (2.5 cm) and 3 ml Luer-Lock Syringe 2: Syringe: 10 ml Luer-Lock 1: Fastener: Catheter Clamp 1: Clamp: Catheter 1: GlideThru™ Peel-Away Sheath: 5 Fr x 2-3/4" (7 cm) Radiopaque over 5 Fr Dilator 2: CSR Wrap 1: Towel: 19-1/2" x 22" (50 x 56 cm) 1: Towel: 24" x 36" (61 x 91 cm) 1: Drape: 36" x 41" (91 x 104 cm) with 3" x 5" (8 x 13 cm) fenestration 1: Drape: 60" x 76" (152 x 193 cm) 5: Gauze Pad: 2" x 2" (5 x 5 cm) 10: Gauze Pad: 4" x 4" (10 x 10 cm) | <ul style="list-style-type: none"> 1: Transducer Cover: 4-12/16" x 95" (12 x 244 cm) 1: Dressing: Tegaderm® I.V. 4" x 6" (10 x 15.5 cm) 1: Dressing: GRIP-LOK™ Catheter Fixation Device 1: Dressing: Steri-Strip® 1/2" x 4" (1.2 x 10 cm) 1: Safety Scalpel: #11 1: Tourniquet⁵ 1: Catheter Trimmer 1: SharpsAway® Disposal Cup 1: SharpsAway II® Locking Disposal Cup 1: Surgical Apparel: Gown 1: Surgical Apparel: Cap 1: Surgical Apparel: Mask² 1: Surgical Apparel: Mask with Shield 1: Paper Tape Measure² 1: HemoHopper® Fluid Receptacle 2: Valve: MicroCLAVE® Neutral Displacement Connector 1: Sterile Procedure Sign² 1: Flow Rate Information Card with Injection Log² 1: Patient ID Card² 1: Chart Sticker² 1: Patient Information Booklet² 1: PICC Line Insertion Checklist² |
|--|--|

All components are CE 0086 unless otherwise noted. Fluid path components are non-pyrogenic. Contains no medication.



Lumen	Priming Volume ¹ (mL)	Flow Rate ¹ (mL/hr)	Pump Flow Rate ² (mL/hr)	MAX Pressure Injection Flow Rate ^{3,4} (mL/sec)
Distal (18 Ga)	0.40	630	3500	4
Proximal (18 Ga)	0.40	720	4090	4

¹ Priming volumes are approximate and are done without injection caps.
² Flow rates were determined using room temperature water, 100 cm head height and represent approximate flow capabilities.
³ Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.
⁴ Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.



- NL** Kit met PICC met twee lumina voor hogedrukinjectie, met Blue FlexTip katheter, positioneringsdraad, prik- en snijbeveiliging en voorzorgsmaatregelen voor maximale steriele barrière
- FR** Kit PICC double lumière pour injection sous pression avec cathéter Blue FlexTip, fil guide, accessoires de protection contre les objets piquants, coupants, tranchants et « précautions barrières maximales »
- DE** Zweilumiges Pressure Injectable PICC (Peripherally Inserted Central Catheter) Kit mit Blue FlexTip Katheter, Platzierungsdraht, Nadel- und Klängen-Entsorgungs-behälter mit Verschluss und maximalen Barriere-Vorsichtsmaßnahmen
- EL** Kit PICC δύο αυτών για έγχυση υπό πίεση, που περιλαμβάνει καθετήρα με άκρο Blue FlexTip, σύρμα τοποθέτησης, διατάξεις προστασίας από αιχμηρά αντικείμενα και μέσα διαφύλαξης μείστων στειρών προφυλάξεων
- IT** Kit con PICC a due lumi per infusione a pressione con catetere Blue FlexTip, funzioni di sicurezza per taglianti e precauzioni per massima barriera di sterilità
- PL** Zestaw cewnika PICC nadającego się do wstrzykiwania pod ciśnieniem, dwukanalowego, z końcówką Blue FlexTip, drutem umieszczającym, zabezpieczeniami przed ukłuciem ostrymi narzędziami i środkami ostrożności zapewniającymi maksymalnie sterylność osłonek
- PT** Kit de cateter central de inserção periférica de dois lúmenes para injeção a pressão com cateter com ponta Blue FlexTip, estilete de introdução, funcionalidades de segurança para objetos cortantes e precauções de protecção integral
- ES** Kit de cateterización central de inserción periférica de dos luces para inyección a presión con catéter con punta Blue FlexTip, fador para colocación, características de seguridad frente a objetos punzantes y precauciones de máxima barrera
- SV** Set med tryckinjicerbar perifert införd central kateter med två lumen och Blue FlexTip, placeringstråd, säkerhetsfunktioner för stickskydd och försiktighetsåtgärder för maximal barriär

Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International, Inc.
2400 Bemville Road
Reading, PA 19605 USA



(01) 0 08 01902 08872 7

SZ-24052-101B (4/15)



Arrow® Pressure Injectable Two-Lumen PICC Kit

REF EU-24052-HPMSB

5 Fr | 2 Lumen | 40 cm catheter length | .018 inch dia. spring-wire guide

