

Dimensions:  
Width: 10-5/8" (27 cm)  
Length: 12-1/2" (31.8 cm)

**Contents:**

- 1: ..... Two-Lumen Indwelling Catheter: 5 Fr x 15-3/4" (40 cm) Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: .. Spring-Wire Guide: .018" (0.46 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
- 1: ..... Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 51-3/16" (130 cm) (Straight Soft Tip on One End — Straight Stiff Tip on Other)
- 1: ..... Introducer Needle: Echogenic 21 Ga x 2-3/4" (7.00 cm) TW
- 1: ..... Syringe: 10 ml Luer-Lock
- 1: ..... Fastener: Catheter Clamp
- 1: ..... Clamp: Catheter
- 1: ..... GlideThru™ Peel-Away Sheath: 5 Fr x 2-3/4" (7 cm) Radiopaque over 5 Fr Dilator
- 1: ..... CSR Wrap
- 1: ..... Dressing: GRIP-LOK® Catheter Fixation Device
- 1: ..... Safety Scalpel: #11
- 1: ..... Catheter Trimmer
- 1: ..... SharpsAway II™ Locking Disposal Cup
- 1: ..... Tourniquet<sup>2</sup>
- 1: ..... Paper Tape Measure<sup>2</sup>
- 2: ... Valve: MicroCLAVE®<sup>1</sup> Neutral Displacement Connector
- 1: ..... Flow Rate Information Card with Injection Log<sup>2</sup>
- 1: ..... Patient ID Card<sup>2</sup>
- 1: ..... Chart Sticker<sup>2</sup>
- 1: ..... Patient Information Booklet<sup>2</sup>

<sup>1</sup> A registered trademark of ICU Medical, Inc.  
<sup>2</sup> Packaged externally to the sterile kit.  
<sup>3</sup> A registered trademark of TIDI.

Fluid path components are non-pyrogenic.  
Contains no medication.

All components are CE 0086 unless otherwise noted.



Arrow International, Inc.  
2400 Bemville Road  
Reading, PA 19605 USA

EC REP Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland


**REF** **EU-24052-IR**

**Pressure Injectable  
Two-Lumen PICC Kit  
with Blue FlexTip® Catheter**

ARROW®



- NL** Kit met PICC met twee lumina voor hogedrukinjectie, met Blue FlexTip catheter
- FR** Kit PICC double lumière pour injection sous pression avec cathéter Blue FlexTip
- DE** Zweilumiges Pressure Injectable PICC (Peripherally Inserted Central Catheter) Kit mit Blue FlexTip Katheter
- EL** Kit PICC δύο αυλών για έγχυση υπό πίεση, που περιλαμβάνει καθετήρα με άκρο Blue FlexTip
- IT** Kit con PICC a due lumi per infusione a pressione con catetere Blue FlexTip
- PL** Zestaw cewnika PICC nadajacego sie do wstrzykiwania pod cisnieniem, dwukanałowego, z końcówką Blue FlexTip
- PT** Kit de cateter central de inserção periférica de dois lúmenes para injeção a pressão com cateter com ponta Blue FlexTip
- ES** Kit de cateterización central de inserción periférica de dos luces para inyección a presión con catéter con punta Blue FlexTip
- SV** Set med tryckinjicerbar perifert införd central kateter med två lumen och Blue FlexTip



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (18 Ga)	0.40	630	3500	4
Proximal (18 Ga)	0.40	720	4090	4

\* Priming volumes are approximate and are done without injection caps.  
† Flow rates were determined using room temperature water, 100 cm head height and represent approximate flow capabilities.  
‡ Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.  
\*\* Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.

**LOT**



(01) 0 08 01902 09303 5

SZ-24052-105A (4/15)



**REF** **EU-24052-IR** **5 Fr** **2 Lumen** **40 cm** catheter length **.018** inch dia. spring-wire guide

Pressure Injectable Two-Lumen PICC Kit with Blue FlexTip® Catheter

