

Dimensions:
Width: 11" (28 cm)
Length: 19" (48.3 cm)



Arrow® Pressure Injectable Multi-Lumen PICC Kit

REF EU-24063-HPMSB

Contents:

- | | |
|---|--|
| <ul style="list-style-type: none"> 1: Multi-Lumen Indwelling Catheter: 6 Fr x 15-3/4" (40 cm) Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, T-Port Connector, Extension Line Clamps, Contamination Guard and Placement Wire 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer™ 1: Filter Needle: 18 Ga TW x 1-1/2" (3.81 cm) with Blunt Fill Tip 1: Introducer Needle: Echogenic 21 Ga x 2-3/4" (7.00 cm) IW 1: Eclipse™ Safety Needle: 22 Ga x 1-5/8" (4 cm) 1: Eclipse™ Safety Needle: 25 Ga x 1" (2.5 cm) and 3 ml Luer-Lock Syringe 2: Syringe: 10 ml Luer-Lock 1: Fastener: Catheter Clamp 1: Clamp: Catheter 1: GlideThru™ Peel-Away Sheath: 6 Fr x 2-3/4" (7 cm) Radiopaque over 6 Fr Dilator 2: CSR Wrap 1: Towel: 19-1/2" x 22" (50 x 56 cm) 1: Towel: 24" x 36" (61 x 91 cm) 1: Drape: 36" x 41" (91 x 104 cm) with 3" x 5" (8 x 13 cm) fenestration 1: Drape: 60" x 76" (152 x 193 cm) 5: Gauze Pad: 2" x 2" (5 x 5 cm) | <ul style="list-style-type: none"> 10: Gauze Pad: 4" x 4" (10 x 10 cm) 1: Transducer Cover: 4-12/16" x 95" (12 x 244 cm) 1: Dressing: Tegaderm® I.V. 4" x 6" (10 x 15.5 cm) 1: Dressing: GRIP-Strip® Catheter Fixation Device 1: Dressing: Steri-Strip® 1/2" x 4" (1.2 x 10 cm) 1: Safety Scalpel: #11 1: Catheter Trimmer 1: SharpsAway® Disposal Cup 1: SharpsAway II™ Locking Disposal Cup 1: Tourniquet† 1: Surgical Apparel: Gown 1: Surgical Apparel: Cap 1: Surgical Apparel: Mask‡ 1: Surgical Apparel: Mask with Shield 1: Paper Tape Measure‡ 1: HemoHopper® Fluid Receptacle 3: Valve: MicroCLAVE® Neutral Displacement Connector 1: Sterile Procedure Sign‡ 1: Flow Rate Information Card with Injection Log‡ 1: Patient ID Card‡ 1: Chart Sticker‡ 1: Patient Information Booklet‡ 1: PICC Line Insertion Checklist‡ |
|---|--|

† A registered trademark of 3M Company.
‡ Packaged externally to the sterile kit.
§ A registered trademark of ICU Medical, Inc.
¶ A trademark of Becton Dickinson and Company.
‡ A registered trademark of TIDI.

All components are CE 0086 unless otherwise noted. Fluid path components are non-pyrogenic. Contains no medication.

Lumen	MAX Pressure			
	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)	Injection Flow Rate§ (mL/sec)
Distal (17 Ga)	0.60	1650	7526	6
Medial (19 Ga)	0.40	550	2736	not rated
Proximal (19 Ga)	0.43	550	2760	not rated

* Priming volumes are approximate and are done without injection caps.
† Flow rates were determined using room temperature water, 100 cm head height and represent approximate flow capabilities.
‡ Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.
§ Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 1/32 cm pressure tubing.



- NL** Kit met PICC met meerdere lumina voor hogedrukinjectie, met Blue FlexTip katheter, positioneringsdraad, prik-en snijbeveiliging en voorzorgsmaatregelen voor maximale steriele barrière
- FR** Kit PICC multi-lumière pour injection sous pression avec cathéter Blue FlexTip, fil guide, accessoires de protection contre les objets piquants, coupants, tranchants et « précautions barrières maximales »
- DE** Mehrlumiges Pressure Injectable PICC (Peripherally Inserted Central Catheter) Kit mit Blue FlexTip Katheter, Platzierungsdraht, Nadel- und Klängen-Entsorgungsbehälter mit Verschluss und maximalen Barriere-Vorsichtsmaßnahmen
- EL** Kit PICC πολλαπλών αυλών για έγχυση υπό πίεση, που περιλαμβάνει καθετήρα με άκρο Blue FlexTip, σύρμα τοποθέτησης, διατάξεις προστασίας από αιχμηρά αντικείμενα και μέσα διαφύλαξης μείγτων στειρών προφυλάξεων
- IT** Kit con PICC a lume multiplo per infusione a pressione con catetere Blue FlexTip, funzioni di sicurezza per taglienti e precauzioni per massima barriera di sterilità
- PL** Zestaw cewnika PICC nadającego się do wstrzykiwania pod ciśnieniem, wielokanalowego, z końcówką Blue FlexTip, drutem umieszczającym, zabezpieczeniami przed ukłuciem ostrymi narzędziami i środkami ostrożności zapewnianymi maksymalnie sterylnej osłone
- PT** Kit de cateter central de inserção periférica de lúmenes múltiplos para injeção a pressão com cateter com ponta Blue FlexTip, estilete de introdução, funcionalidades de segurança para objectos cortantes e precauções de protecção integral
- ES** Kit de cateterización central de inserción periférica múltipluz para inyección a presión con catéter con punta Blue FlexTip, fador para colocación, características de seguridad frente a objetos punzantes y precauciones de máxima barrera
- SV** Set med tryckinjicerbar perifert införd central kateter med flera lumen och Blue FlexTip, placeringstråd, säkerhetsfunktioner för stickskydd och försiktighetsåtgärder för maximal barriär

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International, Inc.
2400 Bemville Road
Reading, PA 19605 USA

LOT



(01) 0 08 01902 08879 6

SZ-24063-101B (4/15)



Arrow® Pressure Injectable Multi-Lumen PICC Kit

REF EU-24063-HPMSB

6 Fr | 3 Lumen | 40 cm catheter length | .018 inch dia. spring-wire guide

