

Dimensions:  
Width: 6-3/4" (17.2 cm)  
Length: 11-5/16" (28.7 cm)

Antimicrobial surfaces utilizing chlorhexidine and silver sulfadiazine. ‡

**REF EU-24701-EN**  
**Central Venous Catheterization Set**  
with Blue FlexTip® ARROWg+ard Blue® Catheter  
**ARROW®**

**Contents:**


- 1: .....Indwelling Catheter: 14 Ga x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard® Antimicrobial Surface Treatment<sup>1</sup>, Integral Extension Line, Extension Line Clamp, Integral Suture Wing
- 1: .....Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: .....Introducer Needle: 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: .....Arrow® Raulerson Spring-Wire Introduction Syringe: 5 ml
- 1: .....Pressure Transduction Probe
- 1: .....Fastener: Catheter Clamp
- 1: .....Clamp: Catheter
- 1: .....Tissue Dilator 7 Fr (2.3 mm) x 4" (10.2 cm)

<sup>1</sup> Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

‡ **Contraindications: The ARROWg+ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).**

Fluid path components are non-pyrogenic. Is not made with natural rubber latex.



Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

**EC REP** Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

- BE** Набор для центральной венозной катетеризации с катетером Blue FlexTip ARROWg+ard Blue
- HR** Komplet za centralnu vensku kateterizaciju s plavim katetrom FlexTip ARROWg+ard Blue
- CS** Centrální žilní katetizační set s katetrem Blue FlexTip ARROWg+ard Blue
- DA** Centralvenekateteriseringssæt med Blue FlexTip ARROWg+ard Blue-kateter
- NL** Centraal veneuze katheterisatieset met Blue FlexTip ARROWg+en Blue katheter
- ET** Tsentraalse veenikateetri komplekt koos Blue FlexTip ARROWg+ard Blue kateetriga
- FI** Keskuslaskimokatetrisarja, jossa sininen FlexTip ARROWg+ard Blue -katetri
- FR** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip ARROWg+ard Blue
- DE** Besteck für zentralen Venenkatheter mit Blue-FlexTip-ARROWg+ard-Blue-Katheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας με Καθέτηρα ARROWg+ard Blue με Blue FlexTip
- HU** Centrális vénás katéterezési szett, Blue FlexTip ARROWg+ard Blue katéterrel
- IS** Miðlaegt æðaleggssett með Blue FlexTip ARROWg+ard Blue æðalegg
- IT** Set per cateterismo della vena centrale con catetere Blue FlexTip ARROWg+ard Blue
- LT** Centrinės venos kateterizavimo rinkinys su „Blue FlexTip ARROWg+ard Blue“ kateteriu
- NO** Sentralt venekateteriseringssett med Blue FlexTip ARROWg+ard Blue-kateter
- PL** Zestaw do cewnikowania żył centralnych z cewnikiem Blue FlexTip ARROWg+ard Blue
- PT** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip, ARROWg+ard Blue
- RO** Set de cateterizare centrală venoasă cu vârf flexibil Blue ARROWg+ard cateter albastru
- RU** Комплект для катетеризации главных вен с катетером Blue FlexTip ARROWg+ard Blue
- SR** ARROWg+ard Blue Catheter komplet za centralnu vensku kateterizaciju sa plavim savitljivim FlexTip vrhom
- SK** Centrálna venózna katetizačná súprava s katétrom ARROWg+ard Blue s nástavcom Blue FlexTip
- SL** Komplet za centralno vensko katetrizacijo s katetrom Blue FlexTip ARROWg+ard Blue
- ES** Conjunto para cateterización venosa central con catéter Blue FlexTip ARROWg+ard Blue
- SV** Set för central venkatetrering med Blue FlexTip ARROWg+ard Blue-kateter
- TR** Santral Venöz Kateterizasyon Seti, Blue FlexTip ARROWg+ard Blue Kateter ile
- UK** Набор для катетеризації центральних вен з катетером Blue FlexTip ARROWg+ard Blue

Catheter Diameter: 14 Ga, 2.2 mm

Priming Volume: 0.55 ml

Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

Flow Rate: 4900 ml/hour

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

SE24701-107A (9/13)

**LOT**



(01) 0 08 01902 08234 3



**REF EU-24701-EN** 14<sup>Ga</sup> 20<sup>cm</sup> catheter length .032<sup>inch dia.</sup> spring-wire guide  
Central Venous Catheterization Set with ARROWg+ard Blue® Catheter

Peel to Open