

REF EU-24701-CVT

Central Venous Catheterization Set

with Blue FlexTip® ARROWg+ard Blue® Catheter

ARROW®

Contents:

- 1:Indwelling Catheter: 14 Ga x 20 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard® Antimicrobial Surface Treatment¹, Integral Extension Line, Extension Line Clamp, Integral Suture Wing
- 1:Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1:Introducer Needle: Echogenic 18 Ga x 2-1/2" (6.35 cm) XTW
- 1:Syringe: 5 ml Luer-Slip
- 1:Tissue Dilator: 7 Fr (2.3 mm) x 10.2 cm
- 1:SecondSite™ Adjustable Hub: Fastener
- 1:SecondSite™ Adjustable Hub: Catheter Clamp
- 1:Scalpel: #11
- 1:ECG Extension Cable

¹ Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

Contraindications: The ARROWg+ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

Fluid path components are non-pyrogenic.



Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Teleflex®
www.teleflex.com

- be** Набор для центральной венозной катетеризации с катетером Blue FlexTip ARROWg+ard Blue
- hr** Komplet za centralnu vensku kateterizaciju s plavim kateterom FlexTip ARROWg+ard Blue
- cs** Centrální žilní katetizační set s katetrem Blue FlexTip ARROWg+ard Blue
- da** Centralvenekateteriseringssæt med Blue FlexTip ARROWg+ard Blue kateter
- nl** Centraal veneuze katheterisatieset met Blue FlexTip ARROWg+ard Blue katheter
- et** Tsentraalse veenikateetri komplekt koos Blue FlexTip ARROWg+ard Blue kateetriga
- fi** Keskuslaskimokatetri sarja, jossa sininen FlexTip ARROWg+ard Blue katetri
- fr** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip ARROWg+ard Blue
- de** Besteck für zentralen Venenkatheter mit Blue FlexTip ARROWg+ard Blue Katheter
- el** Σετ Καθητηριασμού Κεντρικής φλέβας με Καθέτηρα ARROWg+ard Blue με Blue FlexTip
- hu** Centrális vénás katéterezési szett, Blue FlexTip ARROWg+ard Blue katéterrel
- is** Miðlægt æðaleggssett með Blue FlexTip ARROWg+ard Blue æðalegg
- it** Set per cateterismo della vena centrale con catetere Blue FlexTip ARROWg+ard Blue
- lt** Centrinės venos kateterizavimo rinkinys su „Blue FlexTip ARROWg+ard Blue“ kateteriu
- no** Sentralt venekateteriseringssett med Blue FlexTip ARROWg+ard Blue kateter
- pl** Zestaw do cewnikowania żył centralnych z cewnikiem Blue FlexTip ARROWg+ard Blue
- pt** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip, ARROWg+ard Blue
- ro** Set de cateterizare centrală venoasă cu vârf flexibil Blue ARROWg+ard cateter albastru
- ru** Комплект для катетеризации главных вен с катетером Blue FlexTip ARROWg+ard Blue
- sr** ARROWg+ard Blue Catheter komplet za centralnu vensku kateterizaciju sa plavim savitljivim FlexTip vrhom
- sk** Centrálna venózna katetizačná súprava s katétrom ARROWg+ard Blue s nástavcom Blue FlexTip
- sl** Komplet za centralno vensko katetizacijo s katetrom Blue FlexTip ARROWg+ard Blue
- es** Conjunto para cateterización venosa central con catéter Blue FlexTip ARROWg+ard Blue
- sv** Set för central venkateterisering med Blue FlexTip ARROWg+ard Blue kateter
- tr** Santral Venöz Kateterizasyon Seti, Blue FlexTip ARROWg+ard Blue Kateter ile
- uk** Набір для катетеризації центральних вен з катетером Blue FlexTip ARROWg+ard Blue

Catheter Diameter: 14 Ga, 2.2 mm

Priming Volume: 0.55 ml
Priming volumes are approximate and are done without the injection cap.

Flow Rate: 4900 ml/hr
Based on 100 cm head height.

LOT



SZ-24701-111A (5/17)

ARROW



REF

EU-24701-CVT

14^{Ga} 20 cm catheter length .032 inch dia. spring-wire guide

Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

Peel to OPEN