

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

REF EU-24703-CVT
Multi-Lumen Central Venous Catheterization Set
with Blue FlexTip® ARROWg^{ard} Blue® Catheter

ARROW®

Contents:

- 1: Multi-Lumen Indwelling Catheter:
7 Fr x 12" (30 cm) Radiopaque Polyurethane with Blue FlexTip®, ARROWg^{ard}® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with EGG Mark
- 1: Introducer Needle: Echogenic 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Tissue Dilator: 8.5 Fr (2.8 mm) x 10.2 cm
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Scalpel: #11
- 1: ECG Extension Cable
- 2: Injection Site Cap

¹ Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

Contraindications: The ARROWg^{ard} Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

Fluid path components are non-pyrogenic. Not made with natural rubber latex.




Store below 25°C (77°F).
Avoid excessive heat above 40°C (104°F).

Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Многопросветный центральный венозный катетер Набор: Катетер Blue с наконечником Blue FlexTip ARROWg^{ard} Blue Catheter
- HR** Komplet za višelumensku centralno vensku kateterizaciju sa Blue FlexTip ARROWg^{ard} Blue kateterom
- CS** Vícelumenový centrální žilní katetizační set s katetrem Blue FlexTip ARROWg^{ard} Blue
- DA** Multilumen centralvenos kateteriseringsset med Blue FlexTip ARROWg^{ard} Blue kateter
- NL** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip ARROWg^{ard} Blue katheter
- ET** Mitme valendikuga tsentraalse veenikateeteri komplekt koos Blue FlexTip ARROWg^{ard} Blue kateetriga
- F** Monilumeninen keskuskaskimokatetri sarja jossa sininen FlexTip ARROWg^{ard} Blue katetri
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip ARROWg^{ard} Blue
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg^{ard} Blue Katheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας με Πολλαπλούς Αυλούς και Καθετήρα ARROWg^{ard} Blue με Blue FlexTip
- HU** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip ARROWg^{ard} Blue katéterrel
- IS** Fjöllumena holæðaleggsett með Blue FlexTip ARROWg^{ard} Blue æðalegg
- IT** Set per cateterismo della vena centrale multilume con cateter Blue FlexTip ARROWg^{ard} Blue
- LT** Daugiakanaliscentrinės venos kateterizavimo mini rinkinys su „Blue FlexTip ARROWg^{ard} Blue“ antimikrobinu kateteriu
- NO** Multilumen Sentralt venekateteriseringssett med Blue FlexTip ARROWg^{ard} Blue kateter
- PL** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanałowym z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg^{ard} Blue
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip, ARROWg^{ard} Blue
- RO** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip ARROWg^{ard} Blue
- RU** Набор с многопросветным центральным венозным катетером ARROWg^{ard} Blue с кончиком Blue FlexTip
- SR** Višelumenski centralni venski kateterizacioni Set sa Blue FlexTip ARROWg^{ard} Blue kateterom
- SK** Viacľumenová centrálna venózna katetizačná súprava s katétrom Blue FlexTip ARROWg^{ard} Blue
- SL** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip ARROWg^{ard} Blue
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip ARROWg^{ard} Blue
- SV** Set för flerkanalig central venkateterisering med Blue FlexTip ARROWg^{ard} Blue kateter
- TR** ARROWg^{ard} Blue Multi Lumen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Багатопроевітний центральний венозний катетер Набор: катетер Blue с кінчиком Blue FlexTip ARROWg^{ard} Blue Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.55	2400
Medial (18 Ga)	0.45	1070
Proximal (18 Ga)	0.50	1160

* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

5124703-1150 (12/15)

LOT



(01) 0 08 01902 06585 8

ARROW INTERNATIONAL **CVC** **REF EU-24703-CVT** **7 Fr** **3 Lumen** **30 cm** catheter length **.032** inch dia. spring-wire guide **Peel to Open**

Multi-Lumen Central Venous Catheterization Set with ARROWg^{ard} Blue® Catheter