

Dimensions:  
Width: 11" (28 cm)  
Length: 19" (48.3 cm)



## Arrow® Pressure Injectable PICC Kit

REF EU-25041-HPMSB

### Contents:

- |   |   |
|---|---|
| 1: Indwelling Catheter: 4 Fr x 19-5/8" (50 cm)<br>Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, T-Port Connector, Extension Line Clamps, Contamination Guard and Placement Wire | 1: Transducer Cover: 4-12/16" x 95" (12 x 244 cm)             |
| 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer™                                   | 1: Dressing: Tegadem® I.V. 4" x 6" (10 x 15.5 cm)             |
| 1: Filter Needle: 18 Ga TW x 1-1/2" (3.81 cm) with Blunt Fill Tip   | 1: Dressing: GRIP-LOK® Catheter Fixation Device               |
| 1: Introducer Needle: Echogenic 21 Ga x 2-3/4" (7.00 cm) TW   | 1: Dressing: Steri-Strip® 1/2" x 4" (1.2 x 10 cm)             |
| 1: Eclipse™ Safety Needle: 22 Ga x 1-5/8" (4 cm)  | 1: Safety Scalpel: #11  |
| 1: Eclipse™ Safety Needle: 25 Ga x 1" (2.5 cm) and 3 ml Luer-Lock Syringe   | 1: Catheter Trimmer   |
| 2: Syringe: 10 ml Luer-Lock   | 1: SharpsAway® Disposal Cup                                   |
| 1: Fastener: Catheter Clamp   | 1: SharpsAway II™ Locking Disposal Cup                        |
| 1: Clamp: Catheter  | 1: Tourmiquet <sup>2</sup>                                    |
| 1: GlideThru™ Peel-Away Sheath: 4 Fr x 2-3/4" (7 cm) Radiopaque over 4 Fr Dilator   | 1: Surgical Apparel: Gown                                     |
| 2: CSR Wrap   | 1: Surgical Apparel: Cap                                      |
| 1: Towel: 19-1/2" x 22" (50 x 56 cm)  | 1: Surgical Apparel: Mask <sup>2</sup>                        |
| 1: Towel: 24" x 36" (61 x 91 cm)  | 1: Surgical Apparel: Mask with Shield                         |
| 1: Drape: 36" x 41" (91 x 104 cm) with 3" x 5" (8 x 13 cm) fenestration   | 1: Paper Tape Measure <sup>2</sup>                            |
| 1: Drape: 60" x 76" (152 x 193 cm)  | 1: HemoHopper® Fluid Receptacle                               |
| 5: Gauze Pad: 2" x 2" (5 x 5 cm)  | 1: Valve: MicroCLAVE® Neutral Displacement Connector          |
| 10: Gauze Pad: 4" x 4" (10 x 10 cm)   | 1: Sterile Procedure Sign <sup>2</sup>                        |
|   | 1: Flow Rate Information Card with Injection Log <sup>2</sup> |
|   | 1: Patient ID Card <sup>2</sup>                               |
|   | 1: Chart Sticker <sup>2</sup>                                 |
|   | 1: Patient Information Booklet <sup>2</sup>                   |
|   | 1: PICC Line Insertion Checklist <sup>2</sup>                 |

All components are CE 0086 unless otherwise noted. Fluid path components are non-pyrogenic. Contains no medication.

Lumen	Priming Volume* (mL)	Flow Rate <sup>†</sup> (mL/hr)	Pump Flow Rate <sup>‡</sup> (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (18 Ga)	0.50	720	3750	4

\* Priming volumes are approximate and are done without injection caps.  
<sup>†</sup> Flow rates were determined using room temperature water, 100 cm head height and represent approximate flow capabilities.  
<sup>‡</sup> Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.  
 \*\* Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.



- NL** Kit met PICC voor hogedrukinjectie, met Blue FlexTip katheter, positioneringsdraad, prik- en snijbeveiliging en voorzorgsmaatregelen voor maximale steriele barrière
- FR** Kit PICC pour injection sous pression avec cathéter Blue FlexTip, fil guide, accessoires de protection contre les objets piquants, coupants, tranchants et « précautions barrières maximales »
- DE** Pressure Injectable PICC (Peripherally Inserted Central Catheter) Kit mit Blue FlexTip Katheter, Platzierungsdraht, Nadel- und Klagen-Entsorgungsbehälter mit Verschluss und maximalen Barriere-Vorsichtsmaßnahmen
- EL** Kit PICC για έγχυση υπό πίεση, που περιλαμβάνει καθετήρα με άκρο Blue FlexTip, σύρμα τοποθέτησης, διατάξεις προστασίας από αιχμηρά αντικείμενα και μέσα διαφύλαξης μέγιστων στειρών προφυλάξεων
- IT** Kit con PICC per infusione a pressione con catetere Blue FlexTip, funzioni di sicurezza per taglianti e precauzioni per massima barriera di sterilità
- PL** Zestaw nadającego się do wstrzykiwania pod ciśnieniem cewnika PICC z końcówką Blue FlexTip, drutem umieszczającym, zabezpieczeniami przed ukłuciem ostrymi narzędziami i środkami ostrożności zapewniającymi maksymalnie sterylność osłony
- PT** Kit de cateter central de inserção periférica para injeção a pressão com cateter com ponta Blue FlexTip, estilete de introdução, funcionalidades de segurança para objetos cortantes e precauções de protecção integral
- ES** Kit de cateterización central de inserción periférica para inyección a presión con catéter con punta Blue FlexTip, fador para colocación, características de seguridad frente a objetos punzantes y precauciones de máxima barrera
- SV** Set med tryckinjicerbar perifert införd central kateter med Blue FlexTip, placeringstråd, säkerhetsfunktioner för stickskydd och försiktighetsåtgärder för maximal barriär

TELEFLEX Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

LOT



(01) 0 08 01902 08870 3

SZ-25041-101B (5/15)



## Arrow® Pressure Injectable PICC Kit

REF EU-25041-HPMSB

4 Fr

50 cm catheter length

.018 inch dia. spring-wire guide

