

Dimensions:
Width: 10-5/8" (27 cm)
Length: 12-1/2" (31.8 cm)

Contents:

- 1:Indwelling Catheter: 4 Fr x 19-5/8" (50 cm)
Pressure Injectable Radiopaque Polyurethane
with Blue FlexTip[®], Integral Extension Line,
Extension Line Clamp, Integral Suture Wing
- 1: Spring-Wire Guide: .018" (0.46 mm) dia.
x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
- 1: Spring-Wire Guide, Marked: .018" (0.46 mm)
dia. x 51-3/16" (130 cm) (Straight Soft Tip
on One End — Straight Stiff Tip on Other)
- 1: Introducer Needle: Echogenic 21 Ga x 2-3/4"
(7.00 cm) TW
- 1: Syringe: 10 ml Luer-Lock
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: GlideThru™ Peel-Away Sheath: 4 Fr x 2-3/4"
(7 cm) Radiopaque over 4 Fr Dilator
- 1: CSR Wrap
- 1: Dressing: GRIP-LOK^{®3} Catheter Fixation Device
- 1: Safety Scalpel: #11
- 1: Catheter Trimmer
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Tourniquet²
- 1: Paper Tape Measure²
- 1: Valve: MicroCLAVE^{®1} Neutral Displacement Connector
- 1: Flow Rate Information Card with Injection Log²
- 1: Patient ID Card²
- 1: Chart Sticker²
- 1: Patient Information Booklet²

¹ A registered trademark of ICU Medical, Inc.
² Packaged externally to the sterile kit.
³ A registered trademark of TIDI.

Fluid path components are non-pyrogenic.
Contains no medication.

All components are CE 0086 unless otherwise noted.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA


EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

REF
EU-25041-IR

Pressure Injectable PICC Kit with Blue FlexTip[®] Catheter

ARROW[®]

- NL** Kit met PICC voor hogedrukinjectie, met Blue FlexTip katheter
- FR** Kit PICC pour injection sous pression avec cathéter Blue FlexTip
- DE** Pressure Injectable PICC (Peripherally Inserted Central Catheter) Kit mit Blue FlexTip Katheter
- EL** Kit PICC για έγχυση υπό πίεση, που περιλαμβάνει καθετήρα με άκρο Blue FlexTip
- IT** Kit con PICC per infusione a pressione con catetere Blue FlexTip
- PL** Zestaw nadającego się do wstrzykiwania pod ciśnieniem cewnika PICC z końcówką Blue FlexTip
- PT** Kit de cateter central de inserção periférica para injeção a pressão com cateter com ponta Blue FlexTip
- ES** Kit de cateterización central de inserción periférica para inyección a presión con catéter con punta Blue FlexTip
- SV** Set med tryckinjicerbar perifert införd central kateter med Blue FlexTip



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (18 Ga)	0.5	720	3750	4

* Priming volumes are approximate and are done without injection caps.
† Flow rates were determined using room temperature water, 100 cm head height and represent approximate flow capabilities.
‡ Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.
** Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.

LOT



(01) 0 08 01902 09307 3

SZ25041-105A (4/15)



REF

EU-25041-IR

4 Fr

50 cm catheter length

.018 inch dia. spring-wire guide

Pressure Injectable PICC Kit with Blue FlexTip[®] Catheter

