

Dimensions:
Width: 11" (28 cm)
Length: 19" (48.3 cm)



Arrow® Pressure Injectable Multi-Lumen PICC Kit


REF EU-25063-HPMSB

Contents:

- | | |
|--|---|
| 1: Multi-Lumen Indwelling Catheter: 6 Fr x 19.5/8" (50 cm) Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, T-Port Connector, Extension Line Clamps, Contamination Guard and Placement Wire | 1: Transducer Cover: 4-12/16" x 95" (12 x 244 cm) |
| 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer® | 1: Dressing: Tegaderm® I.V. 4" x 6" (10 x 15.5 cm) |
| 1: Filter Needle: 18 Ga TW x 1-1/2" (3.81 cm) with Blunt Fill Tip | 1: Dressing: GRIP-LOK® Catheter Fixation Device |
| 1: Introducer Needle: Echogenic 21 Ga x 2-3/4" (7.00 cm) IW | 1: Dressing: Steri-Strip® 1/2" x 4" (1.2 x 10 cm) |
| 1: Eclipse™ Safety Needle: 22 Ga x 1-5/8" (4 cm) | 1: Safety Scalpel: #11 |
| 1: Eclipse™ Safety Needle: 25 Ga x 1" (2.5 cm) and 3 ml Luer-Lock Syringe | 1: Catheter Trimmer |
| 2: Syringe: 10 ml Luer-Lock | 1: SharpsAway® Disposal Cup |
| 1: Fastener: Catheter Clamp | 1: SharpsAway II™ Locking Disposal Cup |
| 1: Clamp: Catheter | 1: Tourniquet ¹ |
| 1: Glide Thru™ Peel-Away Sheath: 6 Fr x 2-3/4" (7 cm) Radiopaque over 6 Fr Dilator | 1: Surgical Apparel: Gown |
| 2: CSR Wrap | 1: Surgical Apparel: Cap |
| 1: Towel: 19-1/2" x 22" (50 x 56 cm) | 1: Surgical Apparel: Mask ² |
| 1: Towel: 24" x 36" (61 x 91 cm) | 1: Surgical Apparel: Mask with Shield |
| 1: Drape: 36" x 41" (91 x 104 cm) with 3" x 5" (8 x 13 cm) fenestration | 1: Paper Tape Measure ³ |
| 1: Drape: 60" x 76" (152 x 193 cm) | 1: HemolHopper® Fluid Receptacle |
| 5: Gauze Pad: 2" x 2" (5 x 5 cm) | 3: Valve: MicroCLAVE® Neutral Displacement Connector |
| 10: Gauze Pad: 4" x 4" (10 x 10 cm) | 1: Sterile Procedure Sign ² |
| | 1: Flow Rate Information Card with Injection Log ² |
| | 1: Patient ID Card ² |
| | 1: Chart Sticker ² |
| | 1: Patient Information Booklet ² |
| | 1: PICC Line Insertion Checklist ² |

¹ A registered trademark of 3M Company.
² Packaged externally to the sterile kit.
³ A registered trademark of ICU Medical, Inc.
⁴ A trademark of Becton Dickinson and Company.
⁵ A registered trademark of TDI.

All components are CE 0086 unless otherwise noted. Fluid path components are non-pyrogenic. Contains no medication.



Lumen	Priming Volume* (mL)	Flow Rate* (mL/hr)	Pump Flow Rate* (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (17 Ga)	0.65	1320	6765	6
Medial (19 Ga)	0.45	450	2597	not rated
Proximal (19 Ga)	0.46	480	2579	not rated

* Priming volumes are approximate and are done without injection caps.
* Flow rates were determined using room temperature water, 100 cm head height and represent approximate flow capabilities.
* Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.
** Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.



- NL** Kit met PICC met meerdere lumina voor hogedrukinjectie, met Blue FlexTip katheter, positioneringsdraad, prik-en snijbeveiliging en voorzorgsmaatregelen voor maximale steriele barrière
- FR** Kit PICC multi-lumière pour injection sous pression avec cathéter Blue FlexTip, fil guide, accessoires de protection contre les objets piquants, coupants, tranchants et « précautions barrières maximales »
- DE** Mehrlumiges Pressure Injectable PICC (Peripherally Inserted Central Catheter) Kit mit Blue FlexTip Katheter, Platzierungsdraht, Nadel- und Klängen-Entsorgungsbehälter mit Verschluss und maximalen Barriere-Vorsichtsmaßnahmen
- EL** Kit PICC πολλαπλών αυλών για έγχυση υπό πίεση, που περιλαμβάνει καθετήρα με άκρο Blue FlexTip, σύρμα τοποθέτησης, διατάξεις προστασίας από αιχμηρά αντικείμενα και μέσα διαφύλαξης μείστων στειρών προφυλάξεων
- IT** Kit con PICC a lume multiplo per infusione a pressione con catetere Blue FlexTip, funzioni di sicurezza per taglienti e precauzioni per massima barriera di sterilità
- PL** Zestaw cewnika PICC nadającego się do wstrzykiwania pod ciśnieniem, wielokanalowego, z końcówką Blue FlexTip, drutem umieszczającym, zabezpieczeniami przed ukłuciem ostrymi narzędziami i środkami ostrożności zapewnianymi maksymalnie sterylnej osłone
- PT** Kit de cateter central de inserção periférica de lúmenes múltiplos para injeção a pressão com cateter com ponta Blue FlexTip, estilete de introdução, funcionalidades de segurança para objectos cortantes e precauções de protecção integral
- ES** Kit de cateterización central de inserción periférica multiluz para inyección a presión con catéter con punta Blue FlexTip, fador para colocación, características de seguridad frente a objetos punzantes y precauciones de máxima barrera
- SV** Set med tryckinjicerbar perifert införd central kateter med flera lumen och Blue FlexTip, placeringstråd, säkerhetsfunktioner för stickskydd och försiktighetsåtgärder för maximal barriär

ECRIP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International, Inc.
2400 Bemville Road
Reading, PA 19605 USA

LOT



(01) 0 08 01902 08880 2

SZ-25063-101B (5/15)



Arrow® Pressure Injectable Multi-Lumen PICC Kit

REF EU-25063-HPMSB

6 Fr | 3 Lumen | 50 cm catheter length | .018 inch dia. spring-wire guide

