

Dimensions:
Width: 10-5/8" (27 cm)
Length: 12-1/2" (31.8 cm)

Contents:

- 1: Two-Lumen Indwelling Catheter: 5 Fr x 21-21/31" (55 cm) Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, T-Port Connector, Extension Line Clamps, Contamination Guard and Placement Wire
- 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer™
- 1: Introducer Needle: Echogenic 21 Ga x 2-3/4" (7.00 cm) TW
- 1: Syringe: 10 ml Luer-Lock
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: GlideThru™ Peel-Away Sheath: 5 Fr x 2-3/4" (7 cm) Radiopaque over 5 Fr Dilator
- 1: CSR Wrap
- 1: Dressing: GRIP-LOK® Catheter Fixation Device
- 1: Safety Scalpel: #11
- 1: Catheter Trimmer
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Tourniquet²
- 1: Paper Tape Measure²
- 2: ... Valve: MicroCLAVE® Neutral Displacement Connector
- 1: Flow Rate Information Card with Injection Log²
- 1: Patient ID Card²
- 1: Chart Sticker²
- 1: Patient Information Booklet²

¹ A registered trademark of ICU Medical, Inc.
² Packaged externally to the sterile kit.
³ A registered trademark of TIDI

Fluid path components are non-pyrogenic.
Contains no medication.

All components are CE 0086 unless otherwise noted.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

REF **EU-25552-HP**

**Pressure Injectable
Two-Lumen PICC Kit
with Blue FlexTip® Catheter**

ARROW®

- NL** Kit met PICC met twee lumina voor hogedrukinjectie, met Blue FlexTip catheter

FR Kit PICC double lumière pour injection sous pression avec cathéter Blue FlexTip

DE Zweilumiges Pressure Injectable PICC (Peripherally Inserted Central Catheter) Kit mit Blue FlexTip Katheter

EL Kit PICC δύο αυλών για έγχυση υπό πίεση, που περιλαμβάνει καθετήρα με άκρο Blue FlexTip


IT Kit con PICC a due lumi per infusione a pressione con catetere Blue FlexTip

PL Zestaw cewnika PICC nadającego się do wstrzykiwania pod ciśnieniem, dwukanałowego, z końcówką Blue FlexTip

PT Kit de cateter central de inserção periférica de dois lúmenes para injeção a pressão com cateter com ponta Blue FlexTip

ES Kit de cateterización central de inserción periférica de dos lúmenes para inyección a presión con catéter con punta Blue FlexTip

SV Set med tryckinjicerbar perifert införd central kateter med två lumen och Blue FlexTip



LOT

Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (18 Ga)	0.50	480	2934	4
Proximal (18 Ga)	0.50	540	3317	4

* Priming volumes are approximate and are done without injection caps.

† Flow rates were determined using room temperature water, 100 cm head height and represent approximate flow capabilities.

‡ Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.

** Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.



(01) 0 08 01902 09292 2

5225552-103A (3/15)



REF **EU-25552-HP** **5 Fr 2 Lumen 55 cm catheter length .018 inch dia. spring-wire guide**
Pressure Injectable Two-Lumen PICC Kit with Blue FlexTip® Catheter

Radio
Oral