

Dimensions:  
Width: 10-5/8" (27 cm)  
Length: 12-1/2" (31.8 cm)

**Contents:**

- 1: .....Multi-Lumen Indwelling Catheter: 6 Fr x 21-5/8" (55 cm) Pressure Injectable Radiopaque Polyurethane with Blue FlexTip<sup>®</sup>, T-Port Connector, Extension Line Clamps, Contamination Guard and Placement Wire
- 1: .....Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer<sup>™</sup>
- 1: .....Introducer Needle: Echogenic 21 Ga x 2-3/4" (7.00 cm) TW
- 1: ..... Syringe: 10 ml Luer-Lock
- 1: .....Fastener: Catheter Clamp
- 1: .....Clamp: Catheter
- 1: .....GlideThru<sup>™</sup> Peel-Away Sheath: 6 Fr x 2-3/4" (7 cm) Radiopaque over 6 Fr Dilator
- 1: .....CSR Wrap
- 1: .....Dressing: GRIP-LOK<sup>®</sup> Catheter Fixation Device
- 1: .....Safety Scalpel: #11
- 1: .....Catheter Trimmer
- 1: .....SharpsAway II<sup>™</sup> Locking Disposal Cup
- 1: .....Tourniquet<sup>2</sup>
- 1: .....Paper Tape Measure<sup>2</sup>
- 3: .....Valve: MicroCLAVE<sup>®</sup>1 Neutral Displacement Connector
- 1: .....Flow Rate Information Card with Injection Log<sup>2</sup>
- 1: .....Patient ID Card<sup>2</sup>
- 1: .....Chart Sticker<sup>2</sup>
- 1: .....Patient Information Booklet<sup>2</sup>

1 A registered trademark of ICU Medical, Inc.  
2 Packaged externally to the sterile kit.  
3 A registered trademark of TIDI

Fluid path components are non-pyrogenic.  
Contains no medication.

All components are CE 0086 unless otherwise noted.



Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA


EC REP Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

REF
EU-25563-HP

# Pressure Injectable Multi-Lumen PICC Kit with Blue FlexTip<sup>®</sup> Catheter

# ARROW<sup>®</sup>

- NL** Kit met PICC met meerdere lumina voor hogedrukinjectie, met Blue FlexTip katheter
- FR** Kit PICC multi-lumière pour injection sous pression avec cathéter Blue FlexTip
- DE** Viellumiges Pressure Injectable PICC (Peripherally Inserted Central Catheter) Kit mit Blue FlexTip Katheter
- EL** Kit PICC πολλαπλών αυλών για έγχυση υπό πίεση, που περιλαμβάνει καθετήρα με άκρο Blue FlexTip
- IT** Kit con PICC a lume multiplo per infusione a pressione con catetere Blue FlexTip
- PL** Zestaw cewnika PICC nadającego się do wstrzykiwania pod ciśnieniem, wielokanałowego, z końcówką Blue FlexTip
- PT** Kit de cateter central de inserção periférica de lúmenes múltiplos para injeção a pressão com cateter com ponta Blue FlexTip
- ES** Kit de cateterización central de inserción periférica multiluz para inyección a presión con catéter con punta Blue FlexTip
- SV** Set med tryckinjicerbar periferert införd central kateter med flera lumen och Blue FlexTip



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)	MAX Pressure Injection Flow Rate** (mL/sec)
Distal (17 Ga)	0.71	1070	5976	6
Medial (19 Ga)	0.46	360	2246	not rated
Proximal (19 Ga)	0.49	390	2467	not rated

\* Priming volumes are approximate and are done without injection caps.  
† Flow rates were determined using room temperature water, 100 cm head height and represent approximate flow capabilities.  
‡ Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.  
\*\* Pressure injection flow rates are determined at the injector pressure setting of 300 psi maximum using media of 11.8 centipoise viscosity, with 152 cm pressure tubing.

LOT



(01) 0 08 01902 09293 9

SZ25563:103A (3/15)



REF
EU-25563-HP

6 Fr
 3 Lumen
 55 cm catheter length
 .018 inch dia. spring-wire guide

Pressure Injectable Multi-Lumen PICC Kit with Blue FlexTip<sup>®</sup> Catheter

