

Dimensions:  
Width: 6-3/4" (17.2 cm)  
Length: 11-5/16" (28.7 cm)

REF **EU-25802-CVT**

## Two-Lumen Central Venous Catheterization Set

with Blue FlexTip® ARROWg<sup>+</sup>ard Blue® Catheter

# ARROW®

**Contents:**

- 1: .....Two-Lumen Indwelling Catheter:  
8 Fr x 8" (20 cm) Radiopaque  
Polyurethane with Blue FlexTip®,  
ARROWg<sup>+</sup>ard® Antimicrobial Surface  
Treatment<sup>1</sup>, Extension Line Clamps
- 1: .....Spring-Wire Guide, Marked: .032"  
(0.81 mm) dia. x 23-5/8" (60 cm)  
(Straight Soft Tip on One End – "J" Tip  
on Other) with Arrow Advancer™  
with ECG Mark
- 1: .....Introducer Needle: Echogenic  
18 Ga x 2-1/2" (6.35 cm) XTW
- 1: .....Syringe: 5 ml Luer-Slip
- 1: .....Tissue Dilator: 9 Fr (3.0 mm)  
x 10.2 cm
- 1: .....Fastener: Catheter Clamp
- 1: .....Clamp: Catheter
- 1: .....Scalpel: #11
- 1: .....ECG Extension Cable
- 1: .....Injection Site Cap

<sup>1</sup> licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

**Contraindications:** The ARROWg<sup>+</sup>ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

Fluid path components are non-pyrogenic. Not made with natural rubber latex.




Store below 25°C (77°F).  
Avoid excessive heat  
above 40°C (104°F).

Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

**EC REP** Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

- BE** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip ARROWg<sup>+</sup>ard
- HR** Dvolumenski komplet za centralnu vensku kateterizaciju s plavim savitljivim vrhom FlexTip ARROWg<sup>+</sup>ard Blue
- CS** Dvoulumenový centrální žilní katetrizační set s katetrem Blue FlexTip ARROWg<sup>+</sup>ard Blue
- DA** To-lumen centralvenekateteriseringssæt med Blue FlexTip ARROWg<sup>+</sup>ard Blue kateter
- NL** Centraal veneuze dubbellumen katheterisatieset met Blue FlexTip ARROWg<sup>+</sup>ard Blue katheter
- ET** Kahe lumeniga tsentraalveeni kateteriseerimiskomplekt sinise painduva otsaga ARROWg<sup>+</sup>ard Blue kateetriga
- FI** Kaksiluumen-keskuslaskimokatetri setti, jossa Blue FlexTip ARROWg<sup>+</sup>ard Blue katetri
- FR** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip ARROWg<sup>+</sup>ard Blue
- DE** Besteck für zweilumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg<sup>+</sup>ard Blue Katheter
- EL** Σετ Καθητηριασμού Κεντρικής Φλέβας Διπλού Αυλού με Καθέτηρα Blue FlexTip ARROWg<sup>+</sup>ard Blue
- HU** Két lumenű centrális vénás katéterezési készlet Blue FlexTip ARROWg<sup>+</sup>ard Blue katéterrel
- IS** Tveggja-lúmena sett fyrir holæðalegg með Blue FlexTip ARROWg<sup>+</sup>ard Blue æðalegg
- IT** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip ARROWg<sup>+</sup>ard Blue
- LT** Dviejų kanalų centrinės venos kateterizavimo rinkinys su „Blue FlexTip ARROWg<sup>+</sup>ard Blue“ kateteriu
- NO** To-Lumen Sentralt venekateteriseringssett med Blue FlexTip ARROWg<sup>+</sup>ard Blue kateter
- PL** Zestaw do cewnikowania żył centralnych, z 2-kanalowym cewnikiem z końcówką Blue FlexTip i powłoką ARROWg<sup>+</sup>ard Blue
- PT** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip, ARROWg<sup>+</sup>ard Blue
- RO** Trusa de cateterizare venoasă centrală cu două lumene, cu cateter Blue FlexTip ARROWg<sup>+</sup>ard Blue
- RU** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip ARROWg<sup>+</sup>ard
- SR** Dvo-lumenski centralni venski kateterizacijni Set sa Blue FlexTip ARROWg<sup>+</sup>ard Blue kateterom
- SK** Dvojlúmenová katetrizačná súprava na katetrizáciu centrálnej žily s katétrom Blue FlexTip ARROWg<sup>+</sup>ard Blue
- SL** Dvolumenski centralni venski pribor za katetrizacijo s katetrom Blue FlexTip ARROWg<sup>+</sup>ard Blue
- ES** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip ARROWg<sup>+</sup>ard Blue
- SV** Set för tvåkanalyg central venkateterisering med Blue FlexTip ARROWg<sup>+</sup>ard Blue kateter
- TR** Blue FlexTip ARROWg<sup>+</sup>ard Blue Kateterli İki Lümenli Santral Venöz Kateterizasyon Seti
- UK** Двохпросвітний центральный венозный катетер з наконечником Blue FlexTip ARROWg<sup>+</sup>ard



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (14 Ga)	0.75	5600
Proximal (14 Ga)	0.80	5100

\*Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.


LOT



(01) 0 08 01902 06587 2

5725802-1040 (2/15)

ARROW INTERNATIONAL



REF

EU-25802-CVT

8 Fr

2 Lumen

20 cm catheter length

.032 inch dia. spring-wire guide

Two-Lumen Central Venous Catheterization Set with ARROWg<sup>+</sup>ard Blue® Catheter

Peel to Open