

Dimensions:  
Width: 6-3/4" (17.2 cm)  
Length: 11-5/16" (28.7 cm)

REF EU-25853-CVT

## Multi-Lumen Central Venous Catheterization Set

with Blue FlexTip® ARROWg<sup>ard</sup> Blue® Catheter

# ARROW®

**Contents:**

- 1: ..... Multi-Lumen Indwelling Catheter: 8.5 Fr x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, ARROWg<sup>ard</sup>® Antimicrobial Surface Treatment<sup>1</sup>, Extension Line Clamps, Injection Site Caps
- 1: ..... Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1: ..... Introducer Needle: Echogenic 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: ..... Syringe: 5 ml Luer-Slip
- 1: ..... Fastener: Catheter Clamp
- 1: ..... Clamp: Catheter
- 1: ..... Scalpel: #11
- 1: ..... ECG (Electrocardiogram) Extension Cable
- 1: ..... Tissue Dilator

<sup>1</sup> Licensed under US Patent No. 6,706,024.

Fluid path components are non-pyrogenic. Is not made with natural rubber latex.

**Contraindications:** The ARROWg<sup>ard</sup> Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).




Store below 25°C (77°F).  
Avoid excessive heat above 40°C (104°F).

Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

EC REP Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

- BE** Многосветный центральный венозный катетер Набор: Катетер Blue с наконечником Blue FlexTip ARROWg<sup>ard</sup> Blue Catheter
- HR** Komplet za višelumensku centralno vensku kateterizaciju sa Blue FlexTip ARROWg<sup>ard</sup> Blue kateterom
- CS** Višelumenový centrální žilní katetrizační set s katetrem Blue FlexTip ARROWg<sup>ard</sup> Blue
- DA** Multilumen centralvenos kateteriseringsæt med Blue FlexTip ARROWg<sup>ard</sup> Blue kateter
- NL** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip ARROWg<sup>ard</sup> Blue katheter
- ET** Mitme valendikuga tsentraalse veenikateetri komplekt koos Blue FlexTip ARROWg<sup>ard</sup> Blue katetriga
- FI** Monilumeninen keskuskaskimokatetri sarja jossa sininen FlexTip ARROWg<sup>ard</sup> Blue -katetri
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip ARROWg<sup>ard</sup> Blue
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue-FlexTip-ARROWg<sup>ard</sup>-Blue-Katheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας με Πολλαπλούς Ακρούς και Καθετήρα ARROWg<sup>ard</sup> Blue με Blue FlexTip
- HU** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip ARROWg<sup>ard</sup> Blue katéterrel
- IS** Fjöllumena höleðalegsett með Blue FlexTip ARROWg<sup>ard</sup> Blue æðalegg
- IT** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip ARROWg<sup>ard</sup> Blue
- LT** Daugiakanaliscentrinės venos kateterizavimo mini rinkinys su „Blue FlexTip ARROWg<sup>ard</sup> Blue“ antimikrobinū kateteriu
- NO** Multilumen Sentralt venekateteriseringssett med Blue FlexTip ARROWg<sup>ard</sup> Blue-kateter
- PL** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanalowym z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg<sup>ard</sup> Blue
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip, ARROWg<sup>ard</sup> Blue
- RO** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip ARROWg<sup>ard</sup> Blue
- RU** Набор с многосветным центральным венозным катетером ARROWg<sup>ard</sup> Blue с кончиком Blue FlexTip
- SR** Višelumenski centralni venski kateterizacioni Set sa Blue FlexTip ARROWg<sup>ard</sup> Blue kateterom
- SK** Viacľumenová centrálna vložná katetrizačná súprava s katétrom Blue FlexTip ARROWg<sup>ard</sup> Blue
- SL** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip ARROWg<sup>ard</sup> Blue
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip ARROWg<sup>ard</sup> Blue
- SV** Set för flerkanalig central venkateterisering med Blue FlexTip ARROWg<sup>ard</sup> Blue-kateter
- TR** ARROWg<sup>ard</sup> Blue Multi Lumen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Багатопрісвітний центральний венозний катетер Набор: катетер Blue с кінчиком Blue FlexTip ARROWg<sup>ard</sup> Blue Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (16 Ga)	0.47	3640
Medial (14 Ga)	0.84	6640
Proximal (16 Ga)	0.48	3220

\*Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.20 ml.  
†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

5225853-105C (12/13)

LOT



(01) 0 08 01902 06560 5

ARROW  
INTERNATIONAL

REF EU-25853-CVT 8.5<sup>Fr</sup> 3 Lumen 20 cm catheter length .032 inch dia. spring-wire guide  
Multi-Lumen Central Venous Catheterization Set with ARROWg<sup>ard</sup> Blue® Catheter

Peel to Open