

Dimensions:
Width: 6-3/4" (17,2 cm)
Length: 11-5/16" (28,7 cm)

Contents:

- 1: Multi-Lumen Indwelling Catheter:
8.5 Fr x 8" (20 cm) Radiopaque
Polyurethane with Blue FlexTip®,
ARROWg^{ard}® Antimicrobial Surface
Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol,
Marked: .032" (0.81 mm) dia.
x 23-5/8" (60 cm) (Straight Soft Tip
on One End – "J" Tip on Other) with
Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: 18 Ga
x 2-1/2" (6.35 cm) XTW and 5 ml
Luer-Slip Syringe
- 1: Tissue Dilator: 9 Fr (3.0 mm)
x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub:
Catheter Clamp

¹ Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic. Not made with natural rubber latex.

Contraindications: The ARROWg^{ard} Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Store below 25°C (77°F).
Avoid excessive heat
above 40°C (104°F).

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Teleflex®
www.teleflex.com


LOT



REF EU-25853-N
Multi-Lumen Central Venous
Catheterization Set with Blue FlexTip®
ARROWg^{ard} Blue® Catheter

ARROW®

- be** Многопросветный центральный
венозный катетер Набор: Катетер Blue
с наконечником Blue FlexTip ARROWg^{ard}ard
Blue Catheter
- hr** Komplet za višelumensku centralno vensku
kateterizaciju sa Blue FlexTip ARROWg^{ard}ard
Blue kateterom
- cs** Vicelumenový centrální žilní katetrizační set
s katetrem Blue FlexTip ARROWg^{ard}ard Blue
- da** Multilumen centralvenøs kateteriseringssæt
med Blue FlexTip ARROWg^{ard}ard Blue kateter
- nl** Multi-lumen centraal veneuze
katheterisatieset met Blue FlexTip
ARROWg^{ard}ard Blue katheter
- et** Mitme valendikuga tsentraalse veenikateetri
komplekt koos Blue FlexTip ARROWg^{ard}ard
Blue kateetriga
- fi** Monilumeninen keskuskaskimokatetrisarja,
jossa Blue FlexTip ARROWg^{ard}ard Blue katetri
- fr** Jeu de cathétérisme veineux central
à lumières multiples avec cathéter Blue
FlexTip ARROWg^{ard}ard Blue
- de** Besteck für viellumigen, zentralen
Venenkatheter mit Blue FlexTip
ARROWg^{ard}ard Blue Katheter
- el** Σετ Καθετηριασμού Κεντρικής φλέβας
με Πολλαπλούς Αυλούς και Καθετήρα
ARROWg^{ard}ard Blue με Blue FlexTip
- hu** Többlumenű, centrális vénás katéterezési
szett, Blue FlexTip ARROWg^{ard}ard Blue
katéterrel
- is** Fjöl lúmena holæðaleggsett með Blue
FlexTip ARROWg^{ard}ard Blue æðalegg
- it** Set per cateterismo della vena centrale
multilume con catetere Blue FlexTip
ARROWg^{ard}ard Blue
- lt** Daugiakanaliscentrinės venos
kateterizavimo mini rinkinys su „Blue
FlexTip ARROWg^{ard}ard Blue“ antimikrobiniu
kateteriu
- no** Multilumen Sentralt
venekateriseringssett med Blue FlexTip
ARROWg^{ard}ard Blue kateter
- pl** Zestaw do cewnikowania żyły centralnej
z cewnikiem wielokanałowym z końcówką
Blue FlexTip i powłoką antybakteryjną
ARROWg^{ard}ard Blue
- pt** Conjunto de Cateterização Venosa Central
com Lúmen Múltiplo e Cateter Blue
FlexTip, ARROWg^{ard}ard Blue
- ro** Set de cateterizare venoasă centrală
multi-lumen cu cateter Blue FlexTip
ARROWg^{ard}ard Blue
- ru** Набор с многопросветным
центральный венозным катетером
ARROWg^{ard}ard Blue с кончиком
Blue FlexTip
- sr** Višelumenski centralni venski
kateterizacioni Set sa Blue FlexTip
ARROWg^{ard}ard Blue kateterom
- sk** Viacelúmenová centrálna venózná
katetrizačná súprava s katétrom Blue
FlexTip ARROWg^{ard}ard Blue
- sl** Večlumenski centralni venski komplet
za katetrizacijo s katetrom Blue FlexTip
ARROWg^{ard}ard Blue
- es** Conjunto de cateterización venosa central
de luz múltiple con catéter Blue FlexTip
ARROWg^{ard}ard Blue
- sv** Set för flerkanalig central venkatetrisering
med Blue FlexTip ARROWg^{ard}ard Blue
kateter
- tr** ARROWg^{ard}ard Blue Multi Lúmen Santral
Venöz Kateterizasyon Seti, Blue FlexTip
Kateter ile
- uk** Багатопросвітний центральний
венозний катетер Набір: катетер Blue
з кінчиком Blue FlexTip ARROWg^{ard}ard
Blue Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.47	3640
Medial (14 Ga)	0.84	6640
Proximal (16 Ga)	0.48	3220

*Priming volumes are approximate and are done without the injection cap.
†Flow rates are done with normal saline, room temperature,
100 cm head height and represent approximate flow capabilities.

5225853-1088 (11/16)



REF EU-25853-N 8.5 Fr Lumen 20 cm catheter length .032 inch dia. spring-wire guide
Multi-Lumen Central Venous Catheterization Set with ARROWg^{ard}ard Blue® Catheter

Peel to Open