

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

REF **EU-25854-EN**
Quad-Lumen Central Venous
Catheterization Set with Blue FlexTip®
ARROWg+ard Blue® Catheter

ARROW®

Contents:

- 1: Quad-Lumen Indwelling Catheter:
8.5 Fr x 8" (20 cm) Radiopaque
Polyurethane with Blue FlexTip®,
ARROWg+ard® Antimicrobial Surface
Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol,
Marked: .032" (0.81 mm) dia.
x 23-5/8" (60 cm) (Straight Soft Tip
on One End — "J" Tip on Other) with
Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: 18 Ga
x 2-1/2" (6.35 cm) XTW and 5 ml
Arrow® Raulerson Spring-Wire
Introduction Syringe
- 1: Tissue Dilator: 10 Fr (3.5 mm)
x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub:
Catheter Clamp
- 1: Pressure Transduction Probe

¹ Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic. Not made with natural rubber latex.

Contraindications: The ARROWg+ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Store below 25°C (77°F).
Avoid excessive heat
above 40°C (104°F).

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA


EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Teleflex®
www.teleflex.com

- be** Четырехпросветный центральный венозный катетер Набор: Катетер Blue с наконечником Blue FlexTip ARROWg+ard Blue Catheter
- hr** Quad-Lumen Komplet za centralnu vensku kateterizaciju s Blue FlexTip ARROWg+ard Blue kateterom
- cs** Čtyřlumenový centrální žilní katetrizační set s katetrem Blue FlexTip ARROWg+ard Blue
- da** Firdobbelt lumen centralvenekateteriseringsset med Blue FlexTip ARROWg+ard Blue kateter
- nl** Centraal veneuze quad-lumen katheterisatieset met Blue FlexTip ARROWg+ard Blue katheter
- et** Neljavalendikulisea tsentraalse veenikateeteri komplekt koos Blue FlexTip ARROWg+ard Blue kateetriga
- fi** Quad-Lumen keskuslaskimokatetrisarja jossa Blue FlexTip ARROWg+ard Blue katetri
- fr** Jeu de cathétérisme veineux central à quatre lumières avec cathéter Blue FlexTip ARROWg+ard Blue
- de** Besteck für vierlumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg+ard Blue Katheter
- el** Σετ Καθετηριασμού Κεντρικής φλέβας με τετραπλάσιο αυλού και Καθετήρα ARROWg+ard Blue με Blue FlexTip
- hu** Quad-Lumen centrális vénás katérezési szett, Blue FlexTip ARROWg+ard Blue katéterrel
- is** Fjögur lúmena holæðaleggsett með Blue FlexTip ARROWg+ard Blue æðalegg
- it** Set per cateterismo della vena centrale a quattro lumi con catetere Blue FlexTip ARROWg+ard Blue
- lt** Keturių angų centrinės venos kateterizavimo mini rinkinys su „Blue FlexTip ARROWg+ard Blue“ antimikrobiniu kateteriu
- no** Sentralt venekateteriseringssett med fire lumen og Blue FlexTip ARROWg+ard Blue kateter
- pl** Zestaw do cewnikowania żyły centralnej z 4-kanałowym cewnikiem z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg+ard Blue
- pt** Conjunto de Cateterização Venosa Central com Lúmen Quádruplo e Cateter Blue FlexTip, ARROWg+ard Blue
- ro** Set de cateterizare venoasă centrală quad-lumen cu cateter Blue FlexTip ARROWg+ard Blue
- ru** Четырех-просветный набор для центральной венозной катетеризации с катетром Blue FlexTip ARROWg+ard Blue
- sr** Quad-Lumen centralni venški kateterizacijni Set sa Blue FlexTip ARROWg+ard Blue kateterom
- sk** Štvorlúmenová centrálna venózna katetrizačná súprava s katétrom Blue FlexTip ARROWg+ard Blue
- sl** Štiri-lumenski centralni venški komplet za katetrizacijo s katetrom Blue FlexTip ARROWg+ard Blue
- es** Conjunto de cateterización venosa central de cuatro luces con catéter Blue FlexTip ARROWg+ard Blue
- sv** Set för fyrkanalyg central venkateterisering med Blue FlexTip ARROWg+ard Blue kateter
- tr** ARROWg+ard Blue Dörtlü-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- uk** Чотирипросвітний центральний венозний катетер Набор: катетер Blue с кінчиком Blue FlexTip ARROWg+ard Blue Catheter

Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga)	0.43	2700
Medial 1 (14 Ga)	0.54	5100
Medial 2 (18 Ga)	0.36	1400
Proximal (18 Ga)	0.41	1600

* Priming volumes are approximate and are done without the injection cap.
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



REF **EU-25854-EN** 8.5 Fr 4 Lumen 20 cm catheter length .032 inch dia. spring-wire guide
Quad-Lumen Central Venous Catheterization Set with ARROWg+ard Blue® Catheter

5725854-1138 (11/17)

Peel to Open