

Dimensions: 8-5/8" (21.9 cm) x 15-1/2" (39.4 cm)

REF EU-25855-CVT
Multi-Lumen Central Venous Catheterization Set
 with Blue FlexTip® ARROWg⁺ard Blue® Catheter

ARROW®

Contents:

- 1: Multi-Lumen Indwelling Catheter: 8.5 Fr x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, ARROWg⁺ard® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: Echogenic 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Tissue Dilator: 9 Fr (3.0 mm) x 10.2 cm
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Scalpel: #11
- 1: ECG Extension Cable
- 1: Syringe: 5 ml Luer-Slip
- 5: Injection Site Cap

¹ Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

Contraindications: The ARROWg⁺ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

Fluid path components are non-pyrogenic. Not made with natural rubber latex.




Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Arrow International, Inc.
 2400 Bernville Road
 Reading, PA 19605 USA

EC REP Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

- BE** Многопросветный набор для центральной венозной катетеризации с катэтрам Blue FlexTip ARROWg⁺ard Blue Catheter
- HR** Višelumenski Komplet za centralnu vensku kateterizaciju s Blue FlexTip ARROWg⁺ard Blue kateterom
- CS** Vicelumenový centrální žilní katetizační set s katetrem Blue FlexTip ARROWg⁺ard Blue
- DA** Multilumen centralvenekateteriseringsset med Blue FlexTip ARROWg⁺ard Blue kateter
- NL** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip ARROWg⁺ard Blue katheter
- FI** Mitme valendikuga tsentraalse veenikateetri komplekt koost Blue FlexTip ARROWg⁺ard Blue katetriga
- FR** Monilumeninen keskuslaskimokatetrisarja jossa sininen FlexTip ARROWg⁺ard Blue katetri
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip ARROWg⁺ard Blue
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg⁺ard Blue Katheter
- EL** Σετ Καθετηριασμού Κεντρικής φλέβας με Πολλαπλούς Αυλούς και Καθετήρα ARROWg⁺ard Blue με Blue FlexTip
- HU** Többlumenű centrális vénás katéterezési szett, Blue FlexTip ARROWg⁺ard Blue katéterrel
- IS** Fjöllumena holæðaleggsett með Blue FlexTip ARROWg⁺ard Blue æðalegg
- IT** Set per cateterismo della vena centrale a multilume con catetere Blue FlexTip ARROWg⁺ard Blue
- LT** Daugiakanalis centrinės venos kateterizavimo mini rinkinys su „Blue FlexTip ARROWg⁺ard Blue“ antimikrobinu kateteriu
- NO** Multilumen sentralt venekateteriseringssett med og Blue FlexTip ARROWg⁺ard Blue kateter
- PL** Zestaw do cewnikowania żyły centralnej z 5-kanalowym cewnikiem z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg⁺ard Blue
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip, ARROWg⁺ard Blue
- RO** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip ARROWg⁺ard Blue
- RU** Многопросветный центральный венозный катетер с наконечником Blue FlexTip ARROWg⁺ard Blue
- SR** Višelumenski centralni venski kateterizacioni Set sa Blue FlexTip ARROWg⁺ard Blue kateterom
- SK** Vialúmenová centrálna venózna katetrizačná súprava s katetrom Blue FlexTip ARROWg⁺ard Blue
- SL** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip ARROWg⁺ard Blue
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip ARROWg⁺ard Blue
- SV** Set för flerkanalig central venkateterisering med Blue FlexTip ARROWg⁺ard Blue kateter
- TR** Multi Lúmen Santral Venöz Kateterizasyon Seti, Blue FlexTip ARROWg⁺ard Blue Kateter ile
- UK** Багатоприсвітний комплект центрального венозного катетера з катетером Blue с кінчиком Blue FlexTip ARROWg⁺ard Blue Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (16 Ga)	0.63	3240
Medial 1 (14 Ga)	0.78	5720
Medial 2 (18 Ga)	0.38	1070
Medial 3 (18 Ga)	0.40	1220
Proximal (18 Ga)	0.42	1190

*Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.20 ml.

†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

LOT



(01) 0 08 01902 07309 9

SC:25855-1078 (2/15)



REF EU-25855-CVT **8.5^{Fr} 5 Lumen** **20 cm** catheter length **.032** inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg⁺ard Blue® Catheter

Peel to Open