



Two-Lumen Hemodialysis Catheterization Set with Blue FlexTip® ARROWgard Blue® Catheter for High Volume Infusions

REF **CS-26142-F**

14 Fr. | **2** Lumen | **25** cm catheter length | **.038** inch dia. spring-wire guide

Contents:

- 1: Two-Lumen Indwelling Catheter: 14 Fr. x 25 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWgard® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .038" (0.97 mm) dia. x 27-1/2" (70 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 16 Ga. x 1-7/8" (4.8 cm) Radiopaque over 19 Ga. TW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 14 Fr. (4.7 mm) x 14 cm
- 1: Tissue Dilator: 16 Fr. (5.3 mm) x 11.4 cm
- 2: Dust Cap: Non-Vented
- 1: Drape: 24" x 36" (61 cm x 91 cm)
- 1: Scalpel: #11

¹Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.


Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65
Warning: Cancer and Reproductive Harm
www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Contraindications: The Arrow large-bore two-lumen catheter is not designed for long-term hemodialysis or for use in patients with thrombosed vessels.


The ARROWgard Blue antimicrobial surface catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).



LATEX
Not made with natural rubber latex

Store below 25°C (77°F).
Avoid excessive heat above 40°C (104°F).

- fr** Jeu de cathétérisme pour hémodialyse à deux lumières avec cathéter Blue FlexTip ARROWgard Blue pour les perfusions de grand volume
- de** Besteck für zweilumigen Hämodialysekatheter mit Blue-FlexTip-ARROWgard-Blue-Katheter für Infusionen großer Flüssigkeitsmengen
- it** Set per cateterismo in emodialisi a due lumi con catetere Blue FlexTip ARROWgard Blue per volumi di infusione elevati
- pl** Dwukanałowy zestaw do hemodializy z cewnikiem Blue FlexTip ARROWgard Blue do przetaczania dużych objętości płynów infuzyjnych
- pt** Conjunto de Cateterização para Hemodiálise com Lúmen Duplo, Cateter Blue FlexTip, ARROWgard Blue, para Infusões de Grandes Volumes
- ru** Двухпросветный комплект для гемодиализной катетеризации с катетером Blue FlexTip ARROWgard Blue для вливания больших объемов
- sl** Komplet za hemodializno katetrizacijo z dvema svetlinama s katetrom Blue FlexTip ARROWgard Blue za infundiranje velike prostornine
- es** Conjunto de cateterización de dos luces para hemodiálisis con catéter Blue FlexTip ARROWgard Blue para infusiones de gran volumen
- sv** Set för tvåkanalig hemodialyskatetrering med Blue FlexTip ARROWgard Blue-kateter för infusioner med hög volym
- tr** İki Lümenli Hemodiyaliz Kateterizasyon Seti, Blue FlexTip ARROWgard Blue Kateter ile, Yüksek Hacimli İnfüzyonlar için



Lumen	Priming Volume* (mL)	Flow Rate† (mL/min)	Venous Pressure Setting (mmHg)
Distal (10 Ga.)	1.90	250 400	110 180
Proximal (10 Ga.)	1.50	250 400	140 200

* Priming volumes are approximate and are done without accessories.
† Venous Pressure settings are approximate and were determined by varying the blood flow rates from 250 to 400 mL/min in an in vitro experiment using blood analog solution.

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA



teleflex
www.teleflex.com

IS-26142-131A (07/18)



Two-Lumen Hemodialysis Catheterization Set with Blue FlexTip® ARROWgard Blue® Catheter for High Volume Infusions

REF **CS-26142-F**

14 Fr. | **2** Lumen | **25** cm catheter length | **.038** inch dia. spring-wire guide