

Dimensions:  
Width: 6-3/4" (17.2 cm)  
Length: 11-5/16" (28.7 cm)

**REF EU-26553-N**  
**Multi-Lumen Central Venous Catheterization Set**  
with Blue FlexTip® ARROWg<sup>+</sup>ard Blue® Catheter

# ARROW®

**Contents:**

- 1: ..... Multi-Lumen Indwelling Catheter: 5.5 Fr x 5-1/8" (13 cm) Radiopaque Polyurethane with Blue FlexTip®, ARROWg<sup>+</sup>ard® Antimicrobial Surface Treatment<sup>1</sup>, Extension Line Clamps
  - 1: ..... Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
  - 1: ..... Catheter: 22 Ga x 1-3/4" (4.45 cm) Radiopaque over 25 Ga, RW Introducer Needle and 5 ml Luer-Slip Syringe
  - 1: ..... Introducer Needle: 20 Ga x 1-1/2" (4 cm) XTW
  - 1: ..... Tissue Dilator 6 Fr (2.0 mm) x 1-5/8" (4.1 cm)
  - 1: ..... Tissue Dilator 6 Fr (2.0 mm) x 3" (7.6 cm)
  - 1: ..... Fastener: Catheter Clamp
  - 1: ..... Clamp: Catheter
  - 1: ..... Safety Scalpel: #11
  - 1: ... SharpsAway II™ Locking Disposal Cup
- <sup>1</sup> Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

**Contraindications: The ARROWg<sup>+</sup>ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).**

Fluid path components are non-pyrogenic. Is not made with natural rubber latex.




Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

**EC REP** Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

- BE** Многопросветный центральный венозный катетер Набор: Катетер Blue с наконечником Blue FlexTip ARROWg<sup>+</sup>ard Blue Catheter
- HR** Komplet za višelumensku centralno vensku kateterizaciju sa Blue FlexTip ARROWg<sup>+</sup>ard Blue kateterom
- CS** Vícélumenový centrální žilní katetizační set s katetrem Blue FlexTip ARROWg<sup>+</sup>ard Blue
- DA** Multilumen centralvenøs kateteriseringsæt med Blue FlexTip ARROWg<sup>+</sup>ard Blue kateter
- NL** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip ARROWg<sup>+</sup>ard Blue katheter
- ET** Mitme valendikuga tsentraalse veenikateetri komplekt koos Blue FlexTip ARROWg<sup>+</sup>ard Blue kateetriga
- FI** Monilumeninen keskuslaskimokatetrisarja, jossa sininen FlexTip ARROWg<sup>+</sup>ard Blue-katetri
- FR** Jeu de cathétérisme veineux central à lumières multiples avec cathétre Blue FlexTip ARROWg<sup>+</sup>ard Blue
- DE** Besteck für viellumigen, zentralen Venenkatheter mit Blue-FlexTip-ARROWg<sup>+</sup>ard-Blue-Katheter
- EL** Σετ Καθητηριασμού Κεντρικής φλέβας με Πολλαπλούς Αυτούς και Καθητήρα ARROWg<sup>+</sup>ard Blue με Blue FlexTip
- HU** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip ARROWg<sup>+</sup>ard Blue katéterrel
- IS** Fjöllumena holæðaleggsett með Blue FlexTip ARROWg<sup>+</sup>ard Blue æðalegg
- IT** Set per catererismo della vena centrale multilume con catetere Blue FlexTip ARROWg<sup>+</sup>ard Blue
- LT** Daugiakanaliscentrinës venos kateterizavimo mini rinkinys su „Blue FlexTip ARROWg<sup>+</sup>ard Blue“ antimikrobinu kateteriu
- NO** Multilumen Sentralt venekateteriseringssett med Blue FlexTip ARROWg<sup>+</sup>ard Blue-kateter
- PL** Zestaw do cewnikowania żyły centralnej z cewnikami wielokanalowymi z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg<sup>+</sup>ard Blue
- PT** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip, ARROWg<sup>+</sup>ard Blue
- RO** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip ARROWg<sup>+</sup>ard Blue
- RU** Набор с многопросветным центральным венозным катетером ARROWg<sup>+</sup>ard Blue с кончиком Blue FlexTip
- SR** Višelumenski centralni venski kateterizacijski set sa Blue FlexTip ARROWg<sup>+</sup>ard Blue kateterom
- SK** Vialúmenová centrálna venózna katetrizačná súprava s katétrom Blue FlexTip ARROWg<sup>+</sup>ard Blue
- SL** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip ARROWg<sup>+</sup>ard Blue
- ES** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip ARROWg<sup>+</sup>ard Blue
- SV** Set för flerkanalig central venkatetrering med Blue FlexTip ARROWg<sup>+</sup>ard Blue-kateter
- TR** ARROWg<sup>+</sup>ard Blue Multi Lumen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- UK** Багатопроектний центральний венозний катетер Набор: катетер Blue с кілочком Blue FlexTip ARROWg<sup>+</sup>ard Blue Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (20 Ga)	0.30	1200
Medial (22 Ga)	0.31	800
Proximal (22 Ga)	0.34	900

\*Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



(01) 0 08 01902 08761 2

SZ26553102A (1/14)



**REF EU-26553-N** 5.5 Fr 3 Lumen 13 cm catheter length .018 inch dia. spring-wire guide  
Multi-Lumen Central Venous Catheterization Set with ARROWg<sup>+</sup>ard Blue® Catheter

Peel to Open