

Dimensions:  
Width: 6-3/4" (17,2 cm)  
Length: 11-5/16" (28,7 cm)

REF **EU-26702-CVT**

## Two-Lumen Central Venous Catheterization Set

with Blue FlexTip® ARROWg+ard Blue® Catheter

# ARROW®

**Contents:**

- 1: ..... Two-Lumen Indwelling Catheter:  
7 Fr x 8" (20 cm) Radiopaque  
Polyurethane with Blue FlexTip®,  
ARROWg+ard® Antimicrobial Surface  
Treatment<sup>1</sup>, Extension Line Clamps,  
Injection Site Cap
- 1: ..... Spring-Wire Guide, Marked:  
.032" (0.81 mm) dia. x 23-5/8"  
(60 cm) (Straight Soft Tip  
on One End - "J" Tip on Other)  
with Arrow Advancer™
- 1: ..... Introducer Needle: Echogenic  
18 Ga x 2-1/2" (6.35 cm) XTW
- 1: ..... Syringe: 5 ml Luer-Slip
- 1: ..... Fastener: Catheter Clamp
- 1: ..... Clamp: Catheter
- 1: ..... Scalpel: #11
- 1: ..... ECG (Electrocardiogram)  
Extension Cable
- 1: ..... Tissue Dilator

<sup>1</sup> Licensed under US Patent No. 6,706,024.

Fluid path components are non-pyrogenic.  
Is not made with natural rubber latex.

**Contraindications:** The ARROWg+ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).




Store below 25°C (77°F).  
Avoid excessive heat  
above 40°C (104°F).

Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA

EC REP Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

- BE** Двухпросветный центральный венозный катетер з наконечником Blue FlexTip ARROWg+ard
- HR** Dvolumenski komplet za centralnu vensku kateterizaciju s plavim savitljivim vrhom FlexTip ARROWg+ard Blue
- CS** Dvoulumenový centrální žilní katetizační set s katetrem Blue FlexTip ARROWg+ard Blue
- DA** To-lumen centralvenekateteriseringssæt med Blue FlexTip ARROWg+ard Blue kateter
- NL** Centraal veneuze dubbellumen katheterisatieset met Blue FlexTip ARROWg+ard Blue catheter
- ET** Kahe lumeniga tsentraalveeni kateteriseerimiskomplekt sinise painduva otsaga ARROWg+ard Blue katetriga
- FI** Kaksilumenninen keskislaskimokatetri sarja jossa Blue FlexTip ARROWg+ard Blue-katetri
- FR** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip ARROWg+ard Blue
- DE** Besteck für zweilumigen, zentralen Venenkatheter mit Blue-FlexTip-ARROWg+ard-Blue-Katheter
- EL** Σετ Καθετηριασμού Κεντρικής Φλέβας Διπλού Αυλού με Καθετήρα Blue FlexTip ARROWg+ard Blue
- HU** Két lumenű centrális vénás katéterezési készlet Blue FlexTip ARROWg+ard Blue katéterrel
- IS** Tveggja-lúmena sett fyrir holæðalegg með Blue FlexTip ARROWg+ard Blue æðalegg
- IT** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip ARROWg+ard Blue
- LT** Dviejų kanalų centrinės venos kateterizavimo rinkinys su „Blue FlexTip ARROWg+ard Blue“ kateteriu
- NO** To-Lumen Sentralt venekateteriseringssett med Blue FlexTip ARROWg+ard Blue-kateter
- PL** Zestaw do cewnikowania żył centralnych, z 2-kanalowym cewnikiem z końcówką Blue FlexTip i powłoką ARROWg+ard Blue
- PT** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip, ARROWg+ard Blue
- RO** Trusa de cateterizare venoasă centrală cu două lumene, cu cateter Blue FlexTip ARROWg+ard Blue
- RU** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip ARROWg+ard
- SR** Dvo-lumenski centralni venski kateterizacioni Set sa Blue FlexTip ARROWg+ard Blue kateterom
- SK** Dvojlúmenová katetizačná súprava na katetizáciu centrálnej žily s katétrom Blue FlexTip ARROWg+ard Blue
- SL** Dvolumenski centralni venski pribor za katetizacijo s katetrom Blue FlexTip ARROWg+ard Blue
- ES** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip ARROWg+ard Blue
- SV** Set för tvåkanalyg central venekateterisering med Blue FlexTip ARROWg+ard Blue-kateter
- TR** Blue FlexTip ARROWg+ard Blue Kateterli İki Lümenli Santral Venöz Kateterizasyon Seti
- UK** Двохпросвітний центральний венозний катетер з наконечником Blue FlexTip ARROWg+ard



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (16 Ga)	0.47	3700
Proximal (16 Ga)	0.45	3000

\* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

5726702-104C (12/13)

LOT



(01) 0 08 01902 06582 7

**ARROW**  
INTERNATIONAL

REF **EU-26702-CVT**

**7** Fr **2** Lumen **20** cm catheter length **.032** inch dia. spring-wire guide

Two-Lumen Central Venous Catheterization Set with ARROWg+ard Blue® Catheter

Peel to Open