

Dimensions:
Width: 6-3/4" (17,2 cm)
Length: 11-5/16" (28,7 cm)

REF EU-27702-CVT

Two-Lumen Central Venous Catheterization Set

with Blue FlexTip® ARROWg⁺ard Blue® Catheter

ARROW®

Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr x 8" (20 cm) Radiopaque Polyurethane with Blue FlexTip®, ARROWg⁺ard® Antimicrobial Surface Treatment¹, Extension Line Clamps, Injection Site Caps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1: Introducer Needle: Echogenic 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Fastener: Catheter Clamp
- 1: Clamp: Catheter
- 1: Scalpel: #11
- 1: ECG (Electrocardiogram) Extension Cable
- 1: Tissue Dilator

¹ Licensed under US Patent No. 6,706,024.

Fluid path components are non-pyrogenic. Is not made with natural rubber latex.

Contraindications: The ARROWg⁺ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).




Store below 25°C (77°F).
Avoid excessive heat above 40°C (104°F).

Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

- BE** Двухпросветный центральный венозный катетер з наконечником Blue FlexTip ARROWg⁺ard
- HR** Dvolumenski komplet za centralnu vensku kateterizaciju s plavim savitljivim vrhom FlexTip ARROWg⁺ard Blue
- CS** Dvoulumenový centrální žilní katetizační set s katetrem Blue FlexTip ARROWg⁺ard Blue
- DA** To-lumen centralvenekateteriseringsæt med Blue FlexTip ARROWg⁺ard Blue kateter
- NL** Centraal veneuze dubbellumen katheterisatieset met Blue FlexTip ARROWg⁺ard Blue catheter
- ET** Kahe lumeniga tsentraalveeni kateteriseerimiskomplekt sinise painduva otsaga ARROWg⁺ard Blue kateetriga
- FI** Kaksilumeninen keskislaskimokatetri sarja jossa Blue FlexTip ARROWg⁺ard Blue-katetri
- FR** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip ARROWg⁺ard Blue
- DE** Besteck für zweilumigen, zentralen Venenkatheter mit Blue-FlexTip-ARROWg⁺ard-Blue-Katheter
- EL** Σετ Καθετηριασμού Κεντρικής Φλέβας Διπλού Αυλού με Καθετήρα Blue FlexTip ARROWg⁺ard Blue
- HU** Két lumenű centrális vénás katéterezési készlet Blue FlexTip ARROWg⁺ard Blue kateéterrel
- IS** Tveggja-lúmena sett fyrir hölaðalegg með Blue FlexTip ARROWg⁺ard Blue æðalegg
- IT** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip ARROWg⁺ard Blue
- LT** Dvieju kanalu centrines venos kateterizavimo rinkinys su „Blue FlexTip ARROWg⁺ard Blue“ kateteriu
- NO** To-Lumen Sentralt venekateteriseringssett med Blue FlexTip ARROWg⁺ard Blue-kateter
- PL** Zestaw do cewnikowania żył centralnych, z 2-kanalowym cewnikiem z końcówką Blue FlexTip i powłoką ARROWg⁺ard Blue
- PT** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip, ARROWg⁺ard Blue
- RO** Trusa de cateterizare venoasă centrală cu două lumene, cu cateter Blue FlexTip ARROWg⁺ard Blue
- RU** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip ARROWg⁺ard
- SR** Dvo-lumenski centralno venski set za kateterizaciju sa Blue FlexTip ARROWg⁺ard Blue kateterom
- SK** Dvojlúmenová katetizačná súprava na katetrizáciu centrálnych žíl s katétrom Blue FlexTip ARROWg⁺ard Blue
- SL** Dvolumenski centralni venski pribor za kateterizacijo s katetrom Blue FlexTip ARROWg⁺ard Blue
- ES** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip ARROWg⁺ard Blue
- SV** Set för tvåkanalig central venkateterisering med Blue FlexTip ARROWg⁺ard Blue-kateter
- TR** Blue FlexTip ARROWg⁺ard Blue Kateterli İki Lümenli Santral Venöz Kateterizasyon Seti
- UK** Двохпросвітний центральний венозний катетер з наконечником Blue FlexTip ARROWg⁺ard



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)
Distal (14 Ga)	0.57	5000
Proximal (18 Ga)	0.40	1500

*Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

†Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

SE27702-107C (12/13)

LOT



(01) 0 08 01902 06583 4

ARROW
INTERNATIONAL

REF EU-27702-CVT

7 Fr 2 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Two-Lumen Central Venous Catheterization Set with ARROWg⁺ard Blue® Catheter

Peel to Open