

Dimensions:
Width: 6-3/4" (17.2 cm)
Length: 11-5/16" (28.7 cm)

Contents:

- 1: Two-Lumen Indwelling Catheter:
7 Fr x 8" (20 cm) Radiopaque
Polyurethane with Blue FlexTip®,
ARROWg⁺ard Blue® Antimicrobial Surface
Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked:
.032" (0.81 mm) dia.
x 23-5/8" (60 cm) (Straight Soft Tip
on One End – "J" Tip on Other)
with Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: 18 Ga
x 2-1/2" (6.35 cm) XTW
- 1: Arrow® Raulerson Spring-Wire
Introduction Syringe: 5 ml
- 1: Tissue Dilator: 8.5 Fr (2.8 mm)
x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub:
Catheter Clamp
- 1: Pressure Transduction Probe

¹ Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic. Not made with natural rubber latex.

Contraindications: The ARROWg⁺ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Store below 25°C (77°F).
Avoid excessive heat
above 40°C (104°F).

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA


CE REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Teleflex®
www.teleflex.com

REF **EU-27702-EN**
Two-Lumen Central Venous
Catheterization Set with Blue FlexTip®
ARROWg⁺ard Blue® Catheter

ARROW®

- be** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip ARROWg⁺ard
- hr** Dvolumenski komplet za centralnu vensku kateterizaciju s plavim savitljivim vrhom FlexTip ARROWg⁺ard Blue
- cs** Dvoulumenový centrální žilní katetrizační set s katetrem Blue FlexTip ARROWg⁺ard Blue
- da** To-lumen centralvenekateteriseringssæt med Blue FlexTip ARROWg⁺ard Blue kateter
- nl** Centraal veneuze dubbellumen katheterisatieset met Blue FlexTip ARROWg⁺ard Blue katheter
- et** Kahe lumeniga tsentraalveeni kateteriseerimiskomplekt sinise painduva otsaga ARROWg⁺ard Blue kateetriga
- fi** Kaksilumen-keskuslaskimokatetrisetti, jossa Blue FlexTip ARROWg⁺ard Blue katetri
- fr** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip ARROWg⁺ard Blue
- de** Besteck für zweilumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg⁺ard Blue Katheter
- el** Σετ Καθετηριασμού Κεντρικής Φλέβας Διπλού Αυλού με Καθετήρα Blue FlexTip ARROWg⁺ard Blue
- hu** Két lumenű centrális vénás katéterezési készlet Blue FlexTip ARROWg⁺ard Blue katéterrel
- is** Tveggja-lúmena sett fyrir holæðalegg með Blue FlexTip ARROWg⁺ard Blue æðalegg
- it** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip ARROWg⁺ard Blue
- lt** Dviejų kanalų centrinės venos kateterizavimo rinkinys su „Blue FlexTip ARROWg⁺ard Blue“ kateteriu
- no** To-Lumen Sentralt venekateteriseringssett med Blue FlexTip ARROWg⁺ard Blue kateter
- pl** Zestaw do cewnikowania żył centralnych, z 2-kanalowym cewnikiem z końcówką Blue FlexTip i powłoką ARROWg⁺ard Blue
- pt** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip ARROWg⁺ard Blue
- ro** Trusa de cateterizare venoasă centrală cu două lumene, cu cateter Blue FlexTip ARROWg⁺ard Blue
- ru** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip ARROWg⁺ard
- sr** Dvo-lumenski centralni venski set za kateterizaciju sa Blue FlexTip ARROWg⁺ard Blue kateterom
- sk** Dvojlúmenová katetrizačná súprava na katetrizáciu centrálnej žily s katétrom Blue FlexTip ARROWg⁺ard Blue
- sl** Dvolumenski centralni venski pribor za katetrizacijo s katetrom Blue FlexTip ARROWg⁺ard Blue
- es** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip ARROWg⁺ard Blue
- sv** Set för tvåkanalyg central venkateterisering med Blue FlexTip ARROWg⁺ard Blue kateter
- tr** Blue FlexTip ARROWg⁺ard Blue Kateterli İki Lümenli Santral Venöz Kateterizasyon Seti
- uk** Двохпросвітний центральний венозний катетер з наконечником Blue FlexTip ARROWg⁺ard



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (14 Ga)	0.57	5000
Proximal (18 Ga)	0.40	1500

* Priming volumes are approximate and are done without the injection cap.
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

LOT



5427702-111B (11/16)

ARROW



REF EU-27702-EN 7 Fr 2 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Two-Lumen Central Venous Catheterization Set with ARROWg⁺ard Blue® Catheter

Peel to Open