

ARROW®



Arrowgard Blue PLUS® Multi-Lumen CVC Set

REF **CV-45703-E**

7 Fr. | **3** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowgard Blue PLUS® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

¹Licensed under US Patent Nos. 6,706,024, 6,872,195 and 7,329,412.


All components are CE 0086 unless otherwise noted.



Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Fluid path components are non-pyrogenic.

Contraindications: The Arrowgard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



| Lumen | Priming Volume* (mL) | Flow Rate† (mL/hr) |
|-------------------|----------------------|--------------------|
| Distal (16 Ga.) | 0.44 | 3100 |
| Medial (18 Ga.) | 0.39 | 1500 |
| Proximal (18 Ga.) | 0.39 | 1600 |

* Priming volumes are approximate and are done without accessories.
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

- cs** Vícelumenový CVC set Arrowgard Blue PLUS
- nl** Arrowgard Blue PLUS multilumen centraal-veneuze katheterset
- fr** Jeu de cathéters veineux centraux à lumières multiples Arrowgard Blue PLUS
- de** Set mit mehrlumigen Arrowgard Blue PLUS ZVK
- it** Set CVC multilume Arrowgard Blue PLUS
- pl** Wielokanałowy zestaw do cewnikowania żył centralnych z cewnikiem Arrowgard Blue PLUS
- pt** Conjunto de CVC de Lúmen Múltiplo Arrowgard Blue PLUS
- ru** Комплект для катетеризации главных вен с многопросветным катетером Arrowgard Blue PLUS
- sl** Komplet za centralno vensko katetrizacijo z več svetlinami s katetrom Arrowgard Blue PLUS
- es** Conjunto de CVC de luz múltiple Arrowgard Blue PLUS
- sv** Arrowgard Blue PLUS flerlumen CVK-set
- tr** Multi-Lümen Santral Venöz Kateterizasyon Seti, Arrowgard Blue PLUS Kateter ile

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
CE 0086

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IV-45703-134C (09/19)



**Arrowgard Blue PLUS®
Multi-Lumen CVC Set**
REF **CV-45703-E**

7 Fr. | **3** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

ARROW®