

Dimensions:
Width: 9-1/2" (24 cm)
Length: 12-1/2" (31.8 cm)

Contents:

- 1: Quad-Lumen Indwelling Catheter:
8.5 Fr x 8" (20 cm) Pressure Injectable
Radiopaque Polyurethane with Blue FlexTip®,
ARROWg⁺ard Blue PLUS® Antimicrobial
Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked: .032"
(0.81 mm) dia. x 23-5/8" (60 cm)
(Straight Soft Tip on One End – "J" Tip
on Other) with Arrow Advancer™
with ECG Mark
- 1: Introducer Needle: Echogenic 18 Ga
x 2-1/2" (6.35 cm) XTW
- 1: Arrow® Raulerson Spring-Wire Introduction
Syringe: 5 ml
- 1: Tissue Dilator: 10 Fr (3.5 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Safety Scalpel: #11
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Pressure Transduction Probe

¹ Licensed under US Patent Nos. 6,706,024; 6,872,195 and 7,329,412.

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic. Not made with natural rubber latex.

Contraindications: The ARROWg⁺ard Blue® PLUS antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Store below 25°C (77°F).
Avoid excessive heat
above 40°C (104°F).

Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland


REF **EU-45854-CVCPS**

Pressure Injectable
Quad-Lumen Central Venous
Catheterization Kit ARROWg⁺ard Blue PLUS®

MAX
Pressure Injectable

ARROW®

- BE** Набор «ARROWg⁺ard Blue PLUS» для центральной венозной катетеризации четырех полостей органа, які уперськвеецца над націскам
- HR** Komplet za centralnu vensku kateterizaciju s četiri lumena i tlačnim injektiranjem, ARROWg⁺ard Blue PLUS
- CS** Čtyřlumenový centrální žilní kateterizační kit ARROWg⁺ard Blue PLUS, vhodný pro tlakové injekce
- DA** Fire-lumen centralvenekateteriseringsætt til trykinjektion, ARROWg⁺ard Blue PLUS
- NL** Voor injectiepomp geschikte quadri-lumen centraal veneuze katheterisatiekit, ARROWg⁺ard Blue PLUS
- ET** Survega sisestatat nelja valendikuga tsentraalveeni kateteriseerimise komplekt, ARROWg⁺ard Blue PLUS
- FI** Paineinjektoitava neljälumeninen keskuslaskimokatetri sarja, ARROWg⁺ard Blue PLUS
- FR** Kit de cathétérisme veineux central à quatre lumières, injectable sous pression, ARROWg⁺ard Blue PLUS
- DE** Zentralvenenkatheter-Kit mit Vierfach-Lumen zur Druckinjektion, ARROWg⁺ard Blue PLUS
- EL** Κιτ Καθητηριασμού Κεντρικής Φλέβας Τετραπλού Αυλού με Δυνατότητα Εγχύσης με Πίεση, ARROWg⁺ard Blue PLUS
- HU** ARROWg⁺ard Blue PLUS nyomásalápoló befecskendezést alkalmazó négylumenes központi vénás katéterkészlet
- IS** Fjögurra rása miðlaegt holleggjasett fyrir bláæðar sem er dælanlegt með þrýstingi, ARROWg⁺ard Blue PLUS
- IT** Kit per caterismo a quattro lumi per accesso venoso centrale iniettabile a pressione, ARROWg⁺ard Blue PLUS
- LT** Slėginio švirkštimo keturių kanalų centrinės venos kateterizavimo rinkinys, „ARROWg⁺ard Blue PLUS“
- NO** Trykksatt injiserbart fire dobbelt-lumen sentral venekateteriseringssett, ARROWg⁺ard Blue PLUS
- PL** 4-kanalowy zestaw do cewninkowania żył centralnych ARROWg⁺ard Blue PLUS umożliwiający iniekcje pod wysokim ciśnieniem
- PT** Kit de Cateterização Venosa Central de Lúmen Quadruplo Injectável por Pressão ARROWg⁺ard Blue PLUS
- RO** Trusă de cateterizare venoasă centrală cu patru lumene și injectare sub presiune, ARROWg⁺ard Blue PLUS
- RU** Четырехпросветный набор для центральной венозной катетеризации с возможностью введения препарата под давлением, ARROWg⁺ard Blue PLUS
- SR** Komplet sa ubrizgavanjem pritiska i četiri lumena za centralnu vensku kateterizaciju, ARROWg⁺ard Blue PLUS
- SK** Tlaková štvorlúmenová centrálna venózna kateterizačná súprava, ARROWg⁺ard Blue PLUS
- SL** Komplet za centralno vensko kateterizacijo in tlačno injiciranje s štirimi kanali, ARROWg⁺ard Blue PLUS
- ES** Kit de caterismo central venoso de cuatro lúmenes para inyección a presión ARROWg⁺ard Blue PLUS
- SV** Tryckinjicerbara centralvenösa katetersatser med fyra lumen, ARROWg⁺ard Blue PLUS
- TR** Basınçla Enjekte Edilebilen Dört Lümenli Santral Venöz Kateterizasyon Seti, ARROWg⁺ard Blue PLUS
- UK** Чотирипросвітний набір для центральної венозної катетеризації з можливістю введення препарату під тиском, ARROWg⁺ard Blue PLUS



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)	MAX Pressure Injection Flow Rate‡ (ml/sec)
Distal (16 Ga)	0.40	2750	5
Medial 1 (14 Ga)	0.50	5200	10
Medial 2 (18 Ga)	0.36	1300	not rated
Proximal (18 Ga)	0.39	1400	not rated

* Priming volumes are approximate and are done without the injection caps.
† Based on 100 cm head height.
‡ Pressure injection flow rates are determined using media of 11.8 centipoise (cP) viscosity, with 60" administration set tubing.



(01) 0 08 01902 07451 5

SF-45854-108C (10/15)



REF **EU-45854-CVCPS** **8.5 Fr** **4 Lumen** **20 cm** catheter length **.032** inch dia. spring-wire guide
Pressure Injectable Quad-Lumen Central Venous Catheterization Kit, ARROWg⁺ard Blue PLUS®

Peel to Open