

Dimensions:
Width: 9-1/2" (24 cm)
Length: 12-1/2" (31.8 cm)

Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr x 8" (20 cm) Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, ARROWg⁺ard Blue PLUS® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™ with ECG Mark
- 1: Introducer Needle: Echogenic 18 Ga x 2-1/2" (6.35 cm) XTW
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 ml
- 1: Tissue Dilator: 8.5 Fr (2.8 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Safety Scalpel: #11
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Pressure Transduction Probe

¹ Licensed under US Patent Nos. 6,706,024; 6,872,195 and 7,329,412.

All components are CE 0086 unless otherwise noted.

Fluid path components are non-pyrogenic. Not made with natural rubber latex.

Contraindications: The ARROWg⁺ard Blue PLUS antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Store below 25°C (77°F).
Avoid excessive heat above 40°C (104°F).

Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland


REF EU-47702-CVCPS

Pressure Injectable

**Two-Lumen Central Venous
Catheterization Kit, ARROWg⁺ard Blue PLUS®**

ARROW®

- BE** Набор «ARROWg⁺ard Blue PLUS» для центральной венозной катетеризации двух полостей органа, які упрывскваецца пад націскам
- HR** Komplet za centralnu vensku kateterizaciju s dva lumena i tlačnim injektiranjem, ARROWg⁺ard Blue PLUS
- CS** Dvoulumenový centrální žilní katetizační kit ARROWg⁺ard Blue PLUS, vhodný pro tlakové injekce
- DA** Tolumen centralvenekateteriseringsssæt til trykinjektion, ARROWg⁺ard Blue PLUS
- NL** Voor injectiepompe geschikte dubbel-lumen centraal veneuze katheterisatiekit, ARROWg⁺ard Blue PLUS
- ET** Surveiga sisestavat kahe valendikuga tsentraalveeni kateteriseerimise komplekt, ARROWg⁺ard Blue PLUS
- FI** Paineinjektoitava kaksilumeninen keskussaskimokatetri sarja, ARROWg⁺ard Blue PLUS
- FR** Kit de cathétérisme veineux central à deux lumières, injectable sous pression, ARROWg⁺ard Blue PLUS
- DE** Zentralvenenkatheter-Kit mit Doppel-Lumen zur Druckinjektion, ARROWg⁺ard Blue PLUS
- EL** Kit Καθετηριασμού Κεντρικής Φλέβας Διπλού Αυλού με Δυνατότητα Έγχυσης με Πίεση, ARROWg⁺ard Blue PLUS
- HU** ARROWg⁺ard Blue PLUS nyomásalápú befecskendezést alkalmaso kétlumenes központi vénás katéterkészlet
- IS** Tveggja rása miðlægt holleggjasett fyrir bláæðar sem er dælanlegt með þrýstingi, ARROWg⁺ard Blue PLUS
- IT** Kit per cateterismo a doppio lume per accesso venoso centrale iniettabile a pressione, ARROWg⁺ard Blue PLUS
- LT** Slėginio švirkštimo dviejų kanalų centrinės venos kateterizavimo rinkinys, „ARROWg⁺ard Blue PLUS“
- NO** Trykksatt injiserbart to-lumen sentralt venekateteriseringssett, ARROWg⁺ard Blue PLUS
- PL** Dwukanałowy zestaw do cewnikowania żył centralnych ARROWg⁺ard Blue PLUS umożliwiający iniekcje pod wysokim ciśnieniem
- PT** Kit de Cateterização Venoso Central de Lúmen Duplo Injectável por Pressão ARROWg⁺ard Blue PLUS
- RO** Trusă de cateterizare venoasă centrală cu două lumene și injectare sub presiune, ARROWg⁺ard Blue PLUS
- RU** Двухпросветный набор для центральной венозной катетеризации с возможностью введения препарата под давлением, ARROWg⁺ard Blue PLUS
- SR** Komplet sa ubrizgavanjem pritiska i dva lumena za centralnu vensku kateterizaciju, ARROWg⁺ard Blue PLUS
- SK** Tlaková dvoulúmenová centrálna venózna katetizačná súprava, ARROWg⁺ard Blue PLUS
- SL** Komplet za centralno vensko kateterizacijo in tlačno injiciranje z dvema kanaloma, ARROWg⁺ard Blue PLUS
- ES** Kit de cateterismo central venoso de dos lúmenes para inyección a presión ARROWg⁺ard Blue PLUS
- SV** Tryckinjicerbara centralvenösa katetersatsar med två lumen, ARROWg⁺ard Blue PLUS
- TR** Basınçla Enjekte Edilebilen İki Lümenli Santral Venöz Kateterizasyon Seti, ARROWg⁺ard Blue PLUS
- UK** Двохпросвітний набір для центральної венозної катетеризації з можливістю введення препарату під тиском, ARROWg⁺ard Blue PLUS



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hour)	MAX Pressure Injection Flow Rate‡ (ml/sec)
Distal (14 Ga)	0.57	5700	10
Proximal (18 Ga)	0.40	1100	10

* Priming volumes are approximate and are done without injection caps.
† Based on 100 cm head height.
‡ Pressure injection flow rates are determined using media of 11.8 centipoise (cP) viscosity, with 60" administration set tubing.

LOT



SZ47702-108C (10/15)

REF EU-47702-CVCPS

7 Fr 2 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Pressure Injectable Two-Lumen Central Venous Catheterization Kit, ARROWg⁺ard Blue PLUS®

Peel to Open